



## **SBIG Universal Driver/Library**

**Version 4.76 Build 10**

August 13, 2013



**Santa Barbara Instrument Group  
150 Castilian Drive, Suite 101  
Santa Barbara, CA 93117**

**Phone: (805) 308-6985**

**Web: <[www.sbig.com](http://www.sbig.com)> EMail: <[sbig@sbig.com](mailto:sbig@sbig.com)>**

## Table of Contents

<b>1.</b>	<b>Introduction.....</b>	<b>4</b>
<b>2.</b>	<b>Operating System Specific Instructions .....</b>	<b>5</b>
2.1	Microsoft Windows .....	5
2.2	Macintosh OSX.....	5
2.3	Linux .....	6
<b>3.</b>	<b>Driver Interface .....</b>	<b>6</b>
3.1	Driver Related Commands .....	7
3.1.1	Open Driver .....	7
3.1.2	Close Driver.....	7
3.1.3	Open Device .....	7
3.1.4	Close Device .....	8
3.1.5	Get Driver Info.....	8
3.1.6	Get/Set Driver Handle .....	9
3.2	Exposure Related Commands .....	9
3.2.1	Start Exposure, Start Exposure 2 Commands .....	9
3.2.2	End Exposure .....	12
3.2.3	Start Readout.....	12
3.2.4	Readout Line .....	13
3.2.5	Read Subtract Line .....	14
3.2.6	Dump Lines .....	15
3.2.7	End Readout.....	15
3.2.8	Get Line.....	16
3.3	Temperature Related Commands.....	16
3.3.1	Set Temperature Regulation.....	16
3.3.2	Set Temperature Regulation 2.....	18
3.3.3	Query Temperature Status .....	18
3.4	External Control Commands .....	19
3.4.1	Activate Relay .....	20
3.4.2	Pulse Out.....	20
3.4.4	TX Serial Bytes.....	21
3.4.5	Get Serial Status.....	21
3.4.6	AO Tip Tilt .....	21
3.4.7	AO Center .....	21
3.4.8	AO Set Focus.....	22
3.4.9	AO Delay.....	22
3.4.10	CFW.....	22
3.4.11	Motor Focus.....	24
3.5	General Purpose Commands.....	26
3.5.1	Establish Link.....	26
3.5.2	Get CCD Info .....	26
3.5.3	Get Turbo Status.....	29
3.5.4	Query Command Status .....	29
3.5.5	Miscellaneous Control.....	30
3.5.6	Update Clock.....	31
3.5.7	Read Offset.....	31
3.5.8	Read Offset 2.....	31
3.5.9	Get US Timer.....	32
3.5.10	Set/Get IRQ.....	32
3.5.11	Get Link Status.....	32
3.5.12	Get Error String.....	32
3.5.13	Set Driver Control.....	32
3.5.14	Get Driver Control.....	34
3.5.15	USB AD Control .....	34
3.5.16	Query USB.....	35

3.5.17	<i>Query Ethernet</i> .....	36
3.5.18	<i>Get Pentium Cycle Count</i> .....	36
3.5.19	<i>RW USB I2C</i> .....	37
3.5.20	<i>Bit IO</i> .....	37
3.5.21	<i>Column EEPROM</i> .....	38
3.5.22	<i>Customer Options</i> .....	38
3.5.23	<i>Debug Log</i> .....	39
<b>4.</b>	<b>Windows Based Utility Programs</b> .....	<b>41</b>
4.1	SBIGDriverChecker.exe and SBIGDriverChecker64.exe .....	41
4.2	EthSim2.exe.....	41
4.3	SBIGUDRVJournalRx.exe .....	42
4.4	SetClock.exe .....	42
4.5	GetPortD.exe.....	42
<b>5.</b>	<b>Supporting New Cameras and Accessories</b> .....	<b>43</b>
5.1	Supporting the ST-L.....	43
5.2	Supporting the Single Shot Color Cameras.....	44
5.3	Supporting the ST-402 Camera and CFW-402 Filter Wheel .....	44
5.4	Supporting the CFW-10 Color Filter Wheel .....	45
5.5	Supporting the I <sup>2</sup> C based AO Accessories.....	45
5.6	Supporting the I <sup>2</sup> C based CFW-9, CFW-L8 and CFW-L8G Color Filter Wheels.....	45
5.7	Supporting the I <sup>2</sup> C based Motor Focus Accessory .....	45
5.8	Supporting the I <sup>2</sup> C based low-cost A0-8 .....	45
5.9	Supporting the ST-4000XCM Single Shot Color Camera .....	46
5.10	Supporting the STX Cameras .....	46
5.11	Supporting the ST-i Cameras .....	47
5.12	Supporting the STF-8300 Cameras.....	47
5.13	Supporting the STT Cameras .....	47
<b>6.</b>	<b>Revision History</b> .....	<b>48</b>

## 1. Introduction

This document describes the software interface to Santa Barbara Instrument Group's Universal Driver Library (SBIGUDRV). The SBIGUDRV driver supports all of SBIG's Parallel<sup>1</sup>, Ethernet and USB based cameras and accessories. The driver/library is available for Microsoft Windows, Macintosh OSX<sup>2</sup> and Linux. Windows support is for the 32 bit versions of Windows 95<sup>3</sup> through Windows7 and 64 Bit versions of Windows XP through Windows7. The various development kits and supplemental downloads may contain the following files:

**SBIGUDRV.H** - Use this include file with all programs calling the driver. It includes the function prototype and many struct definitions for interfacing to the driver. You need to set the TARGET variable based upon your target environment:

Win 95/98/ME/NT/2000/XP – ENV\_WIN

Macintosh OSX – ENV\_MACOSX

Linux – ENV\_LINUX

One of the things we do with our development is create an LSBIGUDRV.H file where we set TARGET and then include the SBIGUDRV.H file. By including LSBIGUDRV.H instead of SBIGUDRV.H you can replace the SBIGUDRV.H file when it gets revised without having to edit it every time to reset the TARGET.

**SBIGUDRV.LIB** - This is an import library that you link with your 32 bit Windows program. Include this in your "project" file.

**Windows Development Kit** – Contains everything you need to start developing software for the Windows platform including the Windows Tool described below.

**Mac OSX Development Kit** – Contains a USB driver installation package. Running the installer package installs the USB drivers and a Framework that implements the Universal Driver.

**Linux Development Kit** – Contains a gzipped tar archive with the Universal Driver as a shared library and low-level Parallel Port and USB drivers.

**SBIG Driver Checker** – This is an installer utility program for 32 bit versions of a Windows. A complete description is in Section 4.

**SBIG Driver Checker 64** – This is an installer utility program for 64 bit versions of a Windows. A complete description is in Section 4.

**\Tools** – This directory contains some Windows utility programs that will help you get your custom program up and running with the SBIG Universal Driver/Library. See section 3 for a description of each of these tools.

**Visual C++ Sample Program** – Includes a sample Visual Studio 6 project.

**Visual Basic Sample Programs** – Includes some sample Visual Basic programs.

**SBIG C++ Sample Classes** – Includes sample source code for CSBIGCamera and CSBIGImage classes.

---

<sup>1</sup> Parallel cameras are only supported under 32 bit versions of Windows.

<sup>2</sup> Macintosh OSX does not support parallel port based cameras as these ports do not exist in the Mac world. Users with Parallel Port based cameras can use their cameras on Ethernet equipped Macs with OSX by using SBIG's optional Ethernet to Parallel (E2P) Adapter.

<sup>3</sup> There is no support for USB based cameras under Windows 95 or NT.

## 2. Operating System Specific Instructions

This section contains Operating System specific instruction on how to install the drivers and how to build projects that include the SBIG Universal Driver Library.

### 2.1 Microsoft Windows

SBIG provides the SBIG Driver Checker utility for installing the camera drivers under 32 and 64 bit versions of Windows. This process is documented in the “Installing USB.pdf” file in the DOCS directory. It’s critical under Windows that you install the drivers correctly. Please follow that document to the letter. Once the drivers have been installed it’s a good idea to run SBIG’s CCDOps software to make sure you can talk to your camera. If you don’t have CCDOps you can download it from our web site. This will confirm the proper driver installation and insure your camera is working correctly.

Once the drivers are installed and your camera is working follow the instructions below to add the SBIG Universal Driver to your software project:

#### Visual C

- Include the SBIGUDRV.H (or LSBIGUDRV.H) file in your C/C++ files that need to call the driver to get the function prototype, enumeration declarations and struct definitions.
- Add the SBIGUDRV.LIB file to your project so your project will link with the SBIGUDRV.DLL that was installed into the \Windows\System directory

#### Visual BASIC

- Add the SBIGTypes.bas file to your Visual Basic project for the function declaration, enumeration declarations and type definitions.
- Note that Visual BASIC has no unsigned types so pixel data which the driver treats as an unsigned short (values 0 to 65535) but Visual BASIC treats as signed (values –32768 to 32767) will have to be converted to long. Signed values 0 to 32767 convert directly whereas signed values –32768 to –1 convert to 32768 to 65535 respectively. Thus negative signed short pixel values convert to signed long values by “signed long = 65536 – signed short”.

### 2.2 Macintosh OSX

SBIG provides the SBIG Driver Installer package for installing the low-level USB drivers and the SBIG Universal Driver Framework. Double click this file to install the drivers and library. Once the drivers have been installed it’s a good idea to run SBIG’s CCDOps Lite for OSX software to make sure you can talk to your camera. If you don’t have CCDOps Lite you can download it from our web site. This will confirm the proper driver installation and insure your camera is working correctly.

Once the drivers are installed and your camera is working with CCDOps Lite add:

**#include <SBIGUDrv/SBIGUDrv.h>**

to any source code that calls the SBIG Universal Driver Library to bring in the enumerations, struct definitions and function prototype. Also add to the project the following framework file to get the program to link without error:

**/Library/Frameworks/SBIGUDrv.framework**

## 2.3 Linux

Please see the README.txt file included in the Linux Development kit for detailed instructions on how to install the SBIG Universal Driver Library on your Linux system.

## 3. Driver Interface

The software interface to the Universal Driver Library is through single external function that takes a short integer **Command** and pointers to command **Parameter** and command **Results** structs. The driver acts upon the **Command** and **Parameters** struct and fills in the **Results** struct. The memory allocation for these structs is the responsibility of the calling program.

The C prototype for the function is shown in the SBIGUDRV.H file and takes the form:

short **SBIGUnivDrvCommand** (short **Command**, void **\*Parameters**, void **\*Results**)

where **Command** is the command to be executed and **Parameters** and **Results** are pointers to the structs. The function returns an error code indicating whether the camera was able to initiate or complete the command. Note that enumerated types exist in the SBIGUDRV.H file for the **Command** (PAR\_COMMAND) and function return result (PAR\_ERROR) as well as many of the fields in the various **Parameters** and **Results** structs.

The commands supported by the driver are grouped into the following sections discussed individually below:

- Driver Related Commands
- Exposure Related Commands
- Temperature Related Commands
- External Control Commands
- General Purpose Commands

Getting back to the driver, it is written and documented assuming you are programming and proficient in C. As you can see, the function prototype is a C function (not C++), and as you will see the Parameters and Results parameters will end up being pointers to structs using the following data types within the structs:

LOGICAL	- unsigned short (2 bytes) with 0 = FALSE and 1 = TRUE
enum	- Enumerated unsigned short (2 bytes) with an allowed set of values
short	- signed short (2 bytes)
ushort	- unsigned short (2 bytes)
long	- signed long (4 bytes)
ulong	- unsigned long (4 bytes)

Some commands don't require **Parameters** structs and some don't require **Results** structs. In those cases you should pass a NULL pointer to the driver.

The supported commands are discussed in the sections below. For each command the **Parameters** and **Results** structs are shown except in the case where one or both do not exist. The function error return codes for each of the commands will vary from command to command.

For each command that doesn't operate immediately, a command status is maintained internally by the driver, and can be monitored with the **Query Command Status** command. The command status for each of the commands varies from command to command but in general will be from one of the following:

- 0 = Idle
- 1 = Command In Progress

### 3.1 Driver Related Commands

The Commands in this section are used to open and close the driver and get driver related information. At the application level you must Open the Driver (allowing access to the top level of the driver) and then Open the Device (selecting the hardware interface) in order to communicate with the camera.

#### 3.1.1 Open Driver

The Open Driver command is used to initialize the driver and should be your first call to the driver. It takes no Parameters and returns no Results. Just pass NULL pointers to the Parameters and Results arguments of the SBIGUDrvCommand function when you call it.

#### 3.1.2 Close Driver

The Close Driver command is used to close the driver and should be your last call to the driver. There must be one call to Close Driver for each call to Open Driver. This command takes no Parameters and returns no Results. Just pass NULL pointers to the Parameters and Results arguments of the SBIGUDrvCommand function when you call it.

#### 3.1.3 Open Device

The Open Device command is used to load and initialize the low-level driver. You will typically call this second (after Open Driver).

*Parameters Struct:*

```
struct OpenDeviceParams
{
    enum deviceType (see the SBIG_DEVICE_TYPE enum)
        0 = unused
        1,2,3 – LPT 1,2,3
        0x7F00 – USB
        0x7F01 – Ethernet
        0x7F02 – USB1
        0x7F03 – USB2
        0x7F04 – USB3
        0x7F05 – USB4
        0x7F06 – USB5
        0x7F07 – USB6
        0x7F08 – USB7
        0x7F09 – USB8
    ushort lptBaseAddress – for LPT1,2,3 base port address of the LPT port
    ulong ipAddress - for Ethernet the IP address of the camera/accessory
}
```

This command returns no Results. Just pass a NULL pointer to the Results argument of the SBIGUDrvCommand function when you call it.

*Notes:*

- The **lptBaseAddress** is required for LPT1,2,3 under Windows 95/98/Me. This is typically 0x378 for LPT1 and 0x278 for LPT2 but can vary from machine to machine and can be found from the Device Manager control panel. Under Windows NT/2000/XP you can leave this set to 0 as the driver gets this information from the OS.
- The **ipAddress** is required for Ethernet. Use the four bytes of the long with the most significant byte specifying the first part of the address. For example if the desired IP address is 192.168.0.1 use 0xC0A80001 (0xC0=192, 0xA8 = 168, etc.)
- When using the **Open Device** command for opening USB devices, specifying USB (DEV\_USB = 0x7F00) opens the next available USB device which may not be what you want in situations with multiple USB cameras. In these cases use the **Query USB** command first to see what model cameras are available then specify USB1, USB2, USB3 ... USB8 to open those specific cameras.

**3.1.4 Close Device**

The Close Driver command is used to close the low-level driver. You will typically call this second to last (right before Close Driver). There must be one call for Close Device for every call to Open Device.

The Close Device command takes no Parameters and returns no Results. Just pass NULL pointers to the Parameters and Results arguments of the SBIGUDrvCommand function when you call it.

**3.1.5 Get Driver Info**

The Get Driver Info command is used to determine the version and capabilities of the DLL/Driver. For future expandability this command allows you to request several types of information. Initially the standard request and extended requests will be supported but as the driver evolves additional requests will be added.

*Parameters Struct:*

```
struct GetDriverInfoParams
{
    enum request - type of driver information desired (see the DRIVER_REQUEST enum)
        0 = Standard request
        1 = Extended request
        2 = USB loader request
        3, etc. - reserved for future expansion
}
```

*Standard, Extended and USB Loader Results Struct:*

```
struct GetDriverInfoResults0
{
    ushort version - driver version in BCD with the format XX.XX
    char name[64] - driver name, null terminated string
    ushort maxRequest - maximum request response available from this driver
}
```

The Standard request returns the version and name information for the high level SBIG Universal Driver Library. The Extended request returns the version and name information for the low level LPT or USB driver. With linked USB cameras the USB Loader request returns the version and name information for the USB Loader driver.



### 3.1.6 Get/Set Driver Handle

The Get/Set Driver Handle commands are for use by applications that wish to talk to multiple cameras on various ports at the same time. If your software only wants to talk to one camera at a time you can ignore these commands.

The Get Driver Handle command takes a NULL Parameters pointer and a pointer to a GetDriverHandleResults struct for Results. The Set Driver Handle command takes a pointer to a SetDriverHandleParams struct for Parameters and a NULL pointer for Results. To establish links to multiple cameras do the following sequence:

- Call Open Driver for Camera 1
- Call Open Device for Camera 1
- Call Establish Link for Camera 1
- Call Get Driver Handle and save the result as Handle1
- Call Set Driver Handle with INVALID\_HANDLE\_VALUE in the handle parameter
- Call Open Driver for Camera 2
- Call Open Device for Camera 2
- Call Establish Link for Camera 2
- Call Get Driver Handle and save the result as Handle2

Then, when you want to talk to Camera 1, call Set Driver Handle with Handle1 and when you want to talk to Camera 2, call Set Driver Handle with Handle2. To shut down you must call Set Driver Handle, Close Device and Close Driver in that sequence for each camera.

Each time you call Set Driver Handle with INVALID\_HANDLE\_VALUE you are allowing access to an additional camera up to a maximum of four cameras. These cameras can be on different LPT ports, multiple USB<sup>4</sup> cameras or at different Ethernet addresses. There is a restriction though due to memory considerations. You can only have a single readout in process at a time for all cameras and CCDs within a camera. Readout begins with the Start Readout or Readout Line commands and ends with the End Readout command. If you try to do multiple interleaved readouts the data from the multiple cameras will be commingled. To avoid this, simply readout one camera/CCD at a time in an atomic process.

## 3.2 Exposure Related Commands

The commands in the section are used to initiate, complete or cancel an exposure in the camera. For each exposure the camera needs to be instructed to start the exposure, stop the exposure, and readout the image on a row-by-row basis.

### 3.2.1 Start Exposure, Start Exposure 2 Commands

The Start Exposure and Start Exposure 2 commands are used to initiate an exposure. The application specifies the exposure time, etc. and then monitors the exposure's progress with the Query Command Status command discussed below.

The Start Exposure command is the original command from the days when the cameras did not have image buffers and is now deprecated and may eventually be removed. For the interim it can be used with all cameras before the advent of the STX and STT. The Start Exposure 2 replaces it and can be used with all cameras. It should be used in all new and revised code. In addition this command specifies the readout mode and coordinates.

---

<sup>4</sup> At this time supporting multiple USB cameras simultaneously is possible but not very convenient for the user. The first time you open a USB device you will get the first camera connected to the computer, etc. At some point you would like to be able to enumerate the USB cameras. Use the Query USB Command to do that.

*Start Exposure Parameters Struct:*

```

struct StartExposureParams
{
    enum ccd - the CCD to use in the exposure (see the CCD_REQUEST enum)
        0 = Imaging CCD
        1 = Tracking CCD
        2 = External Tracking CCD in ST-L, STX only
    ulong exposureTime - integration time in hundredths of a second in the least significant 24 bits.
        The most significant 8 bits are bit-flags that modify the exposure as described below.
    enum abgState - antiblooming gate state during integration (see the ABG_STATE7 enum)
        0 = Low during integration (ABG shut off)
        1 = Clocked at low rate, 2 = Clocked at medium rate, 3 = Clocked at high rate
    enum openShutter - (see the SHUTTER_COMMAND enum)
        0=Leave Shutter alone, 1=Open Shutter for Exposure and Close for Readout, 2=Close
        Shutter for Exposure and Readout
}

```

*Start Exposure 2 Parameters Struct:*

```

struct StartExposureParams2
{
    enum ccd - same as StartExposureParams above.
    ulong exposureTime - same as StartExposureParams above.
    enum abgState - same as StartExposureParams above.
    enum openShutter - same as StartExposureParams above.
    enum readoutMode - binning mode for readout (see READOUT_BINNING_MODE enum)
        0 = No binning, high resolution
        1 = 2x2 on-chip binning, medium resolution
        2 = 3x3 on-chip binning, low resolution (ST-7/8/etc/237 only)
        0xNN03 = Nx1 on-chip binning (ST-7/8/etc only)
        0xNN04 = Nx2 on-chip binning (ST-7/8/etc only)
        0xNN05 = Nx3 on-chip binning (ST-7/8/etc only)
        6 = No binning, high resolution (ST-7/8/etc only)
        7 = 2x2 binning with vertical binning off-chip (ST-7/8/etc only)
        8 = 3x3 binning with vertical binning off-chip (ST-7/8/etc only)
        9 = 9x9 binning (ST-7/8/etc only)
        0xNN0A = NxN on-chip binning (STF-8300 only)
    ushort top - top most row to readout (0 based)
    ushort left - left most pixel to readout (0 based)
    ushort height - image height in binned pixels
    ushort width - image width in binned pixels
}

```

The status for this command (from the Query Command Status Command) consists of the following two 2-bit fields:

b<sub>1</sub>b<sub>0</sub> = Imaging CCD Status, 00 - CCD Idle, 10=In Progress, 11=Complete

b<sub>3</sub>b<sub>2</sub> = Internal and External Tracking CCD Status, 00 - CCD Idle, 10=In Progress,  
11=Complete

b<sub>15</sub> = Trigger In Status, 1=Waiting for Trigger In, 0=Not waiting

*Notes:*

- For the ST-7/ST-8/ST-9/ST-10/ST-1K/ST-402 add START\_SKIP\_VDD to the **ccd** parameter to increase the image rep rate. This bypasses the time consuming reduction of the CCD's Vdd which is normally used to reduce the readout amplifier glow for the imaging CCD. You'll get a glow in the upper-left corner of the Imaging CCD but the readout rep rate will be higher. SBIG uses this in the Turbo focus mode.
- For the ST-7/ST-8/ST-9/ST-10/ST-1K the minimum allowable exposure is MIN\_ST7\_EXPOSURE (.12 seconds) when the **openShutter** item is 1. If you ask the driver to make a shorter exposure it will take a .12 second exposure.
- For the ST-402 the minimum allowable exposure is MIN\_402\_EXPOSURE (40 milliseconds) when the **openShutter** item is 1. If you ask the driver to make a shorter exposure it will take a 40 millisecond second exposure.
- With Interline CCD cameras like the ST-2K, STL-11K with their electronic shutter the minimum allowable exposure is 0.01 seconds with Version 4.24 of the library and 0.001 seconds with version 4.27 and later.
- To see if a particular camera supports millisecond resolution exposures use the **GetCCDInfo** command Request 4 and check for and Electronic shutter.
- To get millisecond resolution exposures add EXP\_MS\_EXPOSURE to the **exposureTime** item and set the rest of the item to the exposure time in milliseconds. Millisecond exposures from 1 to 255 are possible. Millisecond exposures longer than 255 are rounded up to the nearest 100<sup>th</sup> of a second and programmed as a standard non-millisecond exposure.
- The maximum exposure is 655.35 seconds for Tracking CCD and 167,777.16 seconds for Imaging CCD.
- The **abgState** only affects the TC211 versions of the Tracking CCD on the ST-7/8/etc. and the Imaging CCD of the PixCel255.
- For the PixCel255, PixCel237, ST-1K and ST-402 you need to specify the Imaging CCD since the camera is not a dual CCD design.
- For the PixCel255 and PixCel237 the **openShutter** parameter is ignored and should be set to 0. Use the Pulse Out command to position the Vane/Filter Wheel.
- Bits b<sub>2</sub>/b<sub>3</sub> of the status indicate the Tracking CCD status. With the ST-L, STX and their two tracking CCDs (internal and external) only one set of status is maintained for ease of backwards compatibility. The moral of the story is always issue Start/Stop exposures in pairs for the Tracking CCDs or else these status bits can get confused.
- The ST-L's (STX's) external shutter in the Remote Tracking Head mimics the internal shutter. When the internal shutter is opened the external shutter is opened and vice versa.
- Adding EXP\_WAIT\_FOR\_TRIGGER\_IN to the **exposureTime** item causes the Universal Driver to wait for a Trigger Input signal from the camera before starting the Exposure. This status of the Trigger In signal is indicated by bit b<sub>15</sub> of the command status.
- When you're using the Start Exposure 2 command make sure you still use the Start Readout command and specify the identical readout parameters in both commands.
- With the STF-8300 when you want to enable the Fast Readout mode which effectively cuts the image download time in half at the expense of slightly increased low-level artifacts add the EXP\_FAST\_READOUT to the **exposureTime** item.
- The ST-i only supports four readout modes: High Res (1x1), Medium Res (2x2) and two By-N modes (1xN, 2xN). There is no Low Res mode (3x3) like on many of the other cameras. This information is provided by the **Get CCD Info** command.

### 3.2.2 End Exposure

The End Exposure command is used after the integration is complete to prepare the CCD for readout or to terminate an exposure prematurely.

*Parameters Struct:*

```
struct EndExposureParams
{
    enum ccd - the CCD to end the exposure (see the CCD_REQUEST enum)
        0 = Imaging CCD
        1 = Tracking CCD
        2 = External Tracking CCD in ST-L, STX only
}
```

*Notes:*

- The End Exposure command must be called at least once for each Start Exposure command issued. Several End Exposure commands can be issued without generating an error.
- For the ST-7/8/etc. the End Exposure command prepares the CCD for readout. This normally involves delaying a period of time waiting for the shutter motor to turn off. You can tell the driver to skip this delay by adding END\_SKIP\_DELAY to the **ccd** enum item in order to increase the image rep rate, but you should do this only when the shutter didn't move for both the light and dark images. This means you issued the Start Exposure command with the **openShutter** item set to 0 (leave shutter alone) for the light image and with the **openShutter** item set to 2 (shutter closed for integration and readout) for the dark frame. This scenario only occurs when you are using the Tracking CCD while the Imaging CCD is integrating.
- With the PixCel255, PixCel237, ST-1K and ST-402 the **ccd** parameter should be set to 0 for the Imaging CCD.
- For compatibility with older cameras and software the STX and STT cameras with its internal frame buffer uses the End Exposure command to merely end the exposure phase of an image and start the readout early. If you want to terminate an STX and/or STT exposure completely either or the **ccd** field with ABORT\_DONT\_END (as defined in the SBIGUDRV.H file) or call the End Readout command after the End Exposure command.

### 3.2.3 Start Readout

The Start Readout command is used to inform the driver about the area you intend to readout in subsequent calls to the Readout Line or Read Subtract Line commands. Calling this command is optional (but suggested) and optimizes the readout throughput for small areas on USB and Ethernet based cameras.

*Parameters Struct:*

```
struct StartReadoutParams
{
    enum ccd - the CCD that will be read out (see the CCD_REQUEST enum)
        0 = Imaging CCD
        1 = Tracking CCD
        2 = External Tracking CCD in ST-L, STX only
    enum readoutMode - binning mode for readout (see READOUT_BINNING_MODE enum)
        0 = No binning, high resolution
        1 = 2x2 on-chip binning, medium resolution
}
```

```

    2 = 3x3 on-chip binning, low resolution (ST-7/8/etc/237 only)
    0xNN03 = Nx1 on-chip binning (ST-7/8/etc only)
    0xNN04 = Nx2 on-chip binning (ST-7/8/etc only)
    0xNN05 = Nx3 on-chip binning (ST-7/8/etc only)
    6 = No binning, high resolution (ST-7/8/etc only)
    7 = 2x2 binning with vertical binning off-chip (ST-7/8/etc only)
    8 = 3x3 binning with vertical binning off-chip (ST-7/8/etc only)
    9 = 9x9 binning (ST-7/8/etc only)
    0xNN0A = NxN on-chip binning (STF-8300 only)
    ushort top – topmost row to readout (0 based)
    ushort left – left most pixel to readout (0 based)
    ushort height – image height
    ushort width – image width
}

```

*Notes:*

- See the notes for the Readout Line command.
- The Start Readout command does not actually readout any pixels. It just tells the driver which pixels you will readout using the Readout Line command.
- The Start Readout command does discard **top** lines from the CCD. You do not need to call Dump Lines prior to calls to Readout Line after using the Start Readout Command.
- Even though you specify **left** and **width** parameters with this command you must pass the same values in the **pixelStart** and **pixelLength** parameters in subsequent calls to the Readout Line command.
- With the PixCel255, PixCel237, ST-1K and ST-402 the **ccd** parameter should be set to 0 for the Imaging CCD.
- The ST-i only supports four readout modes: High Res (1x1), Medium Res (2x2) and two By-N modes (1xN, 2xN). There is no Low Res mode (3x3) like on many of the other cameras. This information is provided by the **Get CCD Info** command.

**3.2.4 Readout Line**

The Readout Line command is used to digitize some or all of the active pixels in a row.

*Parameters Struct:*

```

struct ReadoutLineParams
{
    enum ccd - the CCD to readout (see the CCD_REQUEST enum)
        0 = Imaging CCD
        1 = Tracking CCD
        2 = External Tracking CCD in ST-L, STX only
    enum readoutMode - binning mode for readout (see READOUT_BINNING_MODE enum)
        see Start Readout command above
    ushort pixelStart - left most pixel to readout
    ushort pixelLength - number of pixels to digitize
}

```

*Results Struct:*

Rather than passing a pointer to a Results struct, pass a pointer to the destination array of unsigned short integers where the Readout Line command should place the digitized pixel data.

*Notes:*

- Any arbitrary region can be readout using the Dump Lines and Readout Line commands by varying the **pixelStart** and **pixelLength** parameters.
- On Windows 95/98/Me interrupts are disabled for the duration of the line readout. You may want to use the Update Clock command to resynchronize the system clock after reading out an image.
- The PixCel255 and the TC-211 based Tracking CCDs only support the 1x1 and 2x2 binning modes. The 3x3 binning mode is supported by the Imaging CCD, the TC-237 based Tracking CCD and by the PixCel237 only.
- With the PixCel255, PixCel237, ST-1K and ST-402 the **ccd** parameter needs to be set to 0 for the Imaging CCD.
- When binning modes are used, the **pixelStart** and **pixelLength** parameters are in terms of binned pixels.
- You can get the dimensions of the camera's CCD(s) using the Get CCD Info command.
- The Nx1, Nx2 and Nx3 binning modes of the ST-7/8/etc support variable binning (N=1 thru 255) in the vertical direction. You specify the amount of vertical binning in the most significant byte of the readout mode.
- The 1x1, 2x2 and 3x3 off-chip binning modes offer non-streaked horizontal readout for non-Antiblooming versions of the ST-7/8/etc. These detectors bloom both horizontally and vertically under saturating conditions and these readout modes remove the horizontal blooming. They are not necessary with the antiblooming versions of these cameras and the standard on-chip readout modes can be used.
- Readout mode 9 with 9x9 binning is roughly 3 times faster than 3x3 mode and is intended for a fast Focus/Center mode.
- You can only have a single readout in process at a time for all cameras and CCDs in a camera. Readout begins with the Start Readout or Readout Line commands and ends with the End Readout command. If you try to do multiple interleaved readouts the data from the multiple cameras will be commingled. To avoid this, simply readout one camera/CCD at a time in an atomic process.
- The ST-i only supports four readout modes: High Res (1x1), Medium Res (2x2) and two By-N modes (1xN, 2xN). There is no Low Res mode (3x3) like on many of the other cameras. This information is provided by the **Get CCD Info** command.

**3.2.5 Read Subtract Line**

The Read Subtract Line command is identical to the Readout Line command except that it subtracts the data that is stored in memory prior to the readout from the readout data. The Data stored in the array is:  $\text{Data}[n] = \text{CCD}[n] - \text{Data}[n] + 100$

The subtraction adds a bias of 100 to prevent the data from clipping and makes sure the data doesn't overflow or underflow.

*Parameters Struct:*

```
struct ReadoutLineParams
{
    enum ccd - the CCD to readout (see the CCD_REQUEST enum)
        0 = Imaging CCD
```

```

    1 = Tracking CCD
    2 = External Tracking CCD in ST-L, STX only
enum readoutMode - binning mode for readout (see READOUT_BINNING_MODE enum)
    see Start Readout command above
ushort pixelStart - left most pixel to readout
ushort pixelLength - number of pixels to digitize
}

```

#### Results Struct:

Rather than passing a pointer to a Results struct, pass a pointer to the destination array of unsigned short integers where the Read Subtract Line command should place the digitized pixel data.

#### Notes:

- See the notes for the Readout Line command.
- The data is subtracted in place. The Read Subtract command digitizes a pixel, subtracts the value in the destination array, adds 100 counts to avoid clipping at 0 and then stores that result in the destination array.
- With the PixCel255, PixCel237, ST-1K and ST-402 the **ccd** parameter should be set to 0 for the Imaging CCD.

### 3.2.6 Dump Lines

The Dump Lines command is used to discard all of the active pixels in a row on the CCD. You would use this for example when partial frame readout is desired to discard lines above a desired region.

#### Parameters Struct:

```

struct DumpLinesParams
{
    enum ccd - the CCD to dump lines (see the CCD_REQUEST enum)
        0 = Imaging CCD
        1 = Tracking CCD
        2 = External Tracking CCD in ST-L, STX only
    enum readoutMode - binning mode for readout (see READOUT_BINNING_MODE enum)
        see Start Readout command above
    ushort lineLength - number of lines to dump
}

```

#### Notes:

- See the notes for the Readout Line command.
- Unused rows of pixels can be dumped faster than they can be read out. Using the Dump Lines command for sub-array readout can speed up image throughput.
- With the PixCel255, PixCel237, ST-1K and ST-402 the **ccd** parameter should be set to 0 for the Imaging CCD.

### 3.2.7 End Readout

The End Readout command is used after readout of the CCD is complete to prepare the CCD for the idle state.

#### Parameters Struct:

```

struct EndReadoutParams

```

```
{
    enum ccd - the CCD to end the exposure (see the CCD_REQUEST enum)
        0 = Imaging CCD
        1 = Tracking CCD
        2 = External Tracking CCD in ST-L, STX only
}
```

*Notes:*

- The End Readout command should be called at least once per readout after calls to the Readout Line, Read Subtract Line or Dump Lines command are complete. Several End Readout commands can be issued without generating an error.
- For the ST-7/8/etc the End Readout command prepares the CCD for the idle state. This normally involves turning off the CCD preamp and unfreezing the TE cooler if Auto TE Freeze mode has been enabled (see the Set Temperature Regulation command).
- For the other cameras (PixCel255, PixCel237) the End Readout Command does nothing at the current time. For future compatibility you should call this command at the end of the readout phase.
- With the PixCel255, PixCel237, ST-1K and ST-402 the **ccd** parameter should be set to 0 for the Imaging CCD.

### 3.2.8 Get Line

The current driver does not use this command. It was added in a previous version and never removed. It could be reassigned in the future.

## 3.3 Temperature Related Commands

The commands in this section are used to program or monitor the CCD's temperature regulation. Note that parallel port based cameras contain two temperature-sensing thermistors, one in the housing measuring the ambient temperature and one on the CCD. USB based cameras only have a single thermistor mounted on the CCD. Reading the ambient thermistor on those cameras will return a fixed 25°C

### 3.3.1 Set Temperature Regulation

The Set Temperature Regulation command is used to enable or disable the CCD's temperature regulation. It's an older command and uses funny units for the CCD temperature. The Set Temperature Regulation 2 command described in the next section is easier to use with temperatures stated in Degrees Celsius.

*Parameters Struct:*

```
struct SetTemperatureRegulationParams
{
    enum regulation – (see the TEMPERATURE_REGULATION enum)
        0=regulation off, 1=regulation on, 2=regulation override,
        3=freeze TE cooler, 4=unfreeze TE cooler,
        5=enable auto-freeze, 6=disable auto-freeze
    ushort ccdSetpoint - CCD temperature setpoint in A/D units if regulation on or TE drive level (0-255 = 0-100%) if regulation override
}
```

*Notes:*



- The ccdSetpoint above is in A/D units. To convert from temperature in °C to A/D setpoint units and to convert the thermistor readings from the Query Temperature Status command to °C use the following formulas, noting that the CCD thermistor and Ambient thermistor require different constants:

$$T_0 = 25.0$$

$$R_0 = 3.0$$

$$\text{MAX\_AD} = 4096$$

$$R\_RATIO_{\text{CCD}} = 2.57$$

$$R\_RATIO_{\text{Ambient}} = 7.791$$

$$R\_BRIDGE_{\text{CCD}} = 10.0$$

$$R\_BRIDGE_{\text{Ambient}} = 3.0$$

$$DT_{\text{CCD}} = 25.0$$

$$DT_{\text{Ambient}} = 45.0$$

---

#### Calculation of Setpoint from Temperature T in °C

---

$$r = R_0 \times e^{\left( \frac{\ln(R\_RATIO) \times (T_0 - T)}{DT} \right)}$$

$$\text{setpoint} = \frac{\text{MAX\_AD}}{\left( \frac{R\_BRIDGE}{r} + 1.0 \right)}$$

---

#### Calculation of Temperature T in °C from Setpoint

---

$$r = \frac{R\_BRIDGE}{\left( \frac{\text{MAX\_AD}}{\text{setpoint}} - 1.0 \right)}$$

$$T = T_0 - DT \times \left( \frac{\ln\left(\frac{r}{R_0}\right)}{\ln(R\_RATIO)} \right)$$

- What the heck is the "Freeze" all about?** This only pertains to ST-7/8/etc cameras and freezing the TE cooler means telling the temperature regulation circuitry in the camera to keep the TE cooler power at the same level it's currently at until we come back later and unfreeze it. It's essentially a way of telling the camera to be very "quiet" for a period of time. What you do with the freeze commands is use them to freeze the TE cooler for readout and unfreeze it when you're done with the readout. This insures the absolute lowest noise readout possible. One thing you *don't* want to do is keep the TE frozen for a long period of time if you're not reading out the CCD because then you'll degrade the performance of the temperature regulation circuitry. To take advantage of the freeze feature you first tell the camera to enable temperature regulation with the *regulation* item set to 1 and the *setpoint* item set to the desired setpoint temperature, just like you normally would. This allows the camera to achieve regulation at the setpoint temperature. Then, just before you start the readout you call this command with the *regulation* item set to 3 (the *setpoint* does not matter). After readout call this command with the *regulation* item set to 4 (again, *setpoint* doesn't matter). Don't forget to unfreeze the TE or you'll get poor temperature regulation.

You can also get the camera to automatically freeze and unfreeze the TE cooler for you by calling this command with the *regulation* item set to 5. After doing this the driver will freeze the TE at

the first sign of readout and unfreeze it when you call the End Readout command. If you use the auto-freeze feature don't forget that call to End Readout. Finally to disable the auto-freeze function call this command with the *regulation* item set to 6.

Finally, you can query whether the TE is frozen by logically ANDing the *enabled* item of **Query Temperature Status** command results with the REGULATION\_FROZEN\_MASK from the SBIGUDRV.H header. If it's set the TE is currently frozen (either manually or by the auto-freeze feature).

- The cooling in the Remote Guiding Head mimics the internal cooling. While the Remote Guiding Head has unregulated cooling the cooling is enabled whenever internal TE power is above 0%.

### 3.3.2 Set Temperature Regulation 2

The Set Temperature Regulation 2 command is used to enable or disable the CCD's temperature regulation using temperatures in Degrees C instead of the funny A/D units described above.

*Parameters Struct:*

```
struct SetTemperatureRegulationParams2
{
    enum regulation – (see the TEMPERATURE_REGULATION enum)
        0=regulation off, 1=regulation on, 2=regulation override,
        3=freeze TE cooler, 4=unfreeze TE cooler,
        5=enable auto-freeze, 6=disable auto-freeze
    double ccdSetpoint - CCD temperature setpoint in degrees Celsius.
}
```

*Notes:*

- The ccdSetpoint above is in degrees Celsius. This is much easier than using the **SetTemperatureRegulation** command and A/D units as described above.
- For further ease with temperatures in degrees Celsius see the TEMP\_STATUS\_ADVANCED2 request of the **QueryTemperatureStatus** command below.
- Cameras don't have infinite precision for the CCD Setpoint temperatures. The older cameras like the ST-7 allow roughly 0.5° precision whereas the ST-X & STT allows 1/32 °C precision. This command returns the actual camera setpoint in the ccdSetpoint field.

### 3.3.3 Query Temperature Status

The Query Temperature Status command is used to monitor the CCD's temperature regulation. The original version of this command took no Parameters (a NULL pointer) but the command has been expanded to allow a more user friendly result. If you pass a NULL pointer in the Parameters variable you'll get the classic result. If you pass a pointer to a QueryTemperatureStatusParams struct you'll have access to the expanded results.

*Parameters Struct:*

```
struct QueryTemperatureStatusParams
{
    enum request – (see the TEMP_STATUS_REQUEST enum)
}
```

*Standard Results Struct (request = TEMP\_STATUS\_STANDARD):*

```
struct QueryTemperatureStatusResults
{
    LOGICAL enabled - temperature regulation is enabled when this is TRUE
    ushort ccdSetpoint - CCD temperature or thermistor setpoint in A/D units
    ushort power - this is the power being applied to the TE cooler to maintain temperature
        regulation and is in the range 0 thru 255
    ushort ccdThermistor - this is the CCD thermistor reading in A/D units
    ushort ambientThermistor - this is the ambient thermistor reading in A/D units
}
```

*Standard Results Struct (request = TEMP\_STATUS\_ADVANCED2):*

```
struct QueryTemperatureStatusResults2
{
    LOGICAL coolingEnabled - temperature regulation is enabled when this is TRUE
    ushort fanEnabled - fan state and is one of the following: FS_OFF (off), FS_ON (manual
        control) or FS_AUTOCONTROL (auto speed control)
    double ccdSetpoint - CCD Setpoint temperature in °C
    double imagingCCDTemperature - imaging CCD temperature in degrees °C
    double trackingCCDTemperature - tracking CCD temperature in degrees °C
    double externalTrackingCCDTemperature - external tracking CCD temperature in °C
    double ambientTemperature - ambient camera temperature in °C
    double imagingCCDPower - percent power applied to the imaging CCD TE cooler
    double trackingCCDPower - percent power applied to the tracking CCD TE cooler
    double externalTrackingCCDPower - percent power applied to the external tracking TE cooler
    double heatsinkTemperature - imaging CCD heatsink temperature in °C
    double fanPower - percent power applied to the fan
    double fanSpeed - fan speed in RPM
    double trackingCCDSetpoint - tracking CCD Setpoint temperature in °C
}
```

*Notes:*

- Refer to the Set Temperature Regulation command for the formula to convert between A/D units and degrees C for the thermistor readings in the *QueryTemperatureStatusResults*. The advanced results report temperatures directly in Degrees C and need no conversion.
- You can query whether the TE is frozen (see the Set Temperature Regulation command) by logically anding the *enabled* item of the results with the REGULATION\_FROZEN\_MASK from the SBIGUDRV.H header. If it's set the TE is frozen.
- The **fanSpeed** and **trackingCCDSetpoint** were added later and are only returned with the *TEMP\_STATUS\_ADVANCED2* request. A *TEMP\_STATUS\_ADVANCED* request returns all the other values and doesn't set the fanSpeed and trackingCCDSetpoint elements.

### 3.4 External Control Commands

The commands in this section are used to control the telescope position through the telescope interface or to position the CFW-6A motorized color filter wheel.

### 3.4.1 Activate Relay

The Activate Relay command is used to activate one or more of the telescope control outputs or to cancel an activation in progress.

*Parameters Struct:*

```
struct ActivateRelayParams
{
    ushort tXPlus      - x plus activation duration in hundredths of a second
    ushort tXMinus     - x minus activation duration in hundredths of a second
    ushort tYPlus      - y plus activation duration in hundredths of a second
    ushort tYMinus     - y minus activation duration in hundredths of a second
}
```

The status for this command (from the Query Command Status Command) consists of the following four bit field:

b<sub>3</sub> = +X Relay, 0=Off, 1= Active  
 b<sub>2</sub> = -X Relay, 0=Off, 1= Active  
 b<sub>1</sub> = +Y Relay, 0=Off, 1= Active  
 b<sub>0</sub> = -Y Relay, 0=Off, 1= Active

*Notes:*

- This command can be used to cancel relay activations by setting the appropriate parameters to 0.

### 3.4.2 Pulse Out

The Pulse Out command is used with the ST-7/8/etc to position the CFW-6A/CFW-8 and with the PixCel255 and PixCel237 to position the internal vane/filter wheel.

*Parameters Struct:*

```
struct PulseOutParams
{
    ushort numberPulses - number of pulses to generate (0 thru 255)
    ushort pulseWidth   - width of pulses in units of microseconds with a minimum of 9 microseconds
    ushort pulsePeriod  - period of pulses in units of microseconds with a minimum of 29 plus the
                           pulseWidth microseconds
}
```

The status for this command (from the Query Command Status command) consists of the following bit fields:

b<sub>0</sub> - Normal status, 0 = inactive, 1 = pulse out in progress  
 b<sub>1</sub>-b<sub>3</sub> - PixCel255/237 Filter state, 0=moving, 1-5=at position 1-5, 6=unknown

*Notes:*

- The camera will cease communications while the Pulse Out command is in progress to maintain the best pulse width accuracy. After sending the ACK response the camera will generate the pulses and only when it has finished generating the pulses will it respond to further communications from the PC.

- With the PixCel255/237 for positioning the internal vane/filter wheel you set the numberPulses parameter to a non-zero value (typically 1), the pulseWidth to zero and the pulsePeriod to one of the following values: 0=Leave vane/filter alone, 1-5=Position vane/filter wheel at position 1 thru 5, 6=Stop motor, abort any move in progress, 7=initialize and identify vane/filter wheel.
- On the PixCel255/237 the following filter positions are defined: Position 1 = Clear/Open, Position 2 = Opaque, Position 3 = Red, Position 4 = Green and Position 5 = Blue. Positions 1 and 2 are supported by the vane and positions 1 thru 5 are supported by the filter wheel.
- You find out what type of filter wheel is installed in the PixCel255/237 using the Get CCD Info command with request number 3.
- See the **CFW** command below for a high level API for programming the SBIG color filter wheels.
- The Pulse Out command directed to an ST-L will cause the ST-L's CFW-L to emulate a CFW-8. This makes old software for the CFW-8 work on the CFW-L but new code should use the CFW command to support the CFW-L.

#### 3.4.4 TX Serial Bytes

The TX Serial Bytes command is for internal use by SBIG. It's a very low level version of commands like AO Tip Tilt that are used to send data out the ST-7/8/etc's telescope port to accessories like the AO-7. There's no reason why you should need to use this command. Just use the dedicated commands like AO Tip Tilt.

#### 3.4.5 Get Serial Status

The Get Serial Status command is for internal use by SBIG. It's a very low level version of commands like AO Tip Tilt that are used to send data out the ST-7/8/etc's telescope port to accessories like the AO-7. There's no reason why you should need to use this command. Just use the dedicated commands like AO Tip Tilt.

#### 3.4.6 AO Tip Tilt

The AO Tip Tilt Command is used to position an AO-7 attached to the telescope port of an ST-7/8/etc.

*Parameters Struct:*

```
struct AOTipTiltParams
{
    ushort xDeflection - this is the desired position of the mirror in the X axis
    ushort yDeflection - this is the desired position of the mirror in the Y axis
}
```

*Notes:*

- The range for the X and Y deflection parameters are 0 through 4095. The mirror is centered at 2048, fully to one side at 0 and fully at the other side with 4095.
- While commanding the AOL to position 2048, 2048 will take it to the center of its range, using the AO Center command described below will force the AOL to re-find center using the mechanical home positioning sensors.

#### 3.4.7 AO Center

This command centers the AO attached to the camera. It takes no *Parameters* or *Results* structs and returns an error code. This works with the AO-7 and the AOL but in the case of the AOL this causes the unit to find the center using the mechanical home sensors.

### 3.4.8 AO Set Focus

This command is reserved for future use with motorized focus units. Prototypes of the AO-7 had motorized focus but the feature was removed in the production units. This command is a holdover from that.

### 3.4.9 AO Delay

The AO Delay Command is used to generate millisecond type delays for exposing the Tracking CCD.

*Parameters Struct:*

```
struct AODelayParams
{
    ulong delay - this is the desired delay in microseconds
}
```

*Notes:*

- The computer essentially hangs while waiting for this delay to expire so be careful how you use this command.

### 3.4.10 CFW

The CFW Command is a high-level API for controlling the SBIG color filter wheels. It supports the **CFW-2** (two position shutter wheel in the ST-5C/237), the **CFW-5** (internal color filter wheel for the ST-5C/237), the **CFW-8**, the internal filter wheel (**CFW-L**) in the ST-L Large Format Camera, the internal filter wheel (**CFW-402**) in the ST-402 camera, the old 6-position **CFW-6A**, the 10-position **CFW-10** in both I2C and RS-232 interface modes, the I2C based **CFW-9** and 8-position CFW for the STL (**CFW-L8**), the five (**FW5-STX**) and seven (**FW7-STX**) position CFWs for the STX, the five (**FW5-8300**) and eight (**FW8-8300**) position CFWs for the ST-8300 and the eight (**FW8-STT**) position CFW for the STT cameras.

*Parameters Struct:*

```
struct CFWParams
{
    enum cfwModel – (see the CFW_MODEL_SELECT enum)
        0=Unknown, 1=CFW-2, 2=CFW-5,
        3=CFW-8, 4=CFW-L, 5=CFW-402,
        6=Auto detect, 7=CFW-6A, 8=CFW-10,
        9=RS232 based CFW-10, 10=CFW-9, 11=Standard CFW-L8,
        12=Custom CFW-L8, 13=CFW-1603, 14=FW5-STX,
        15=FW5-8300, 16=FW8-8300, 17=FW7-STX, 18=FW8-STT

    enum cfwCommand – (see the CFW_COMMAND enum)
        0=Query, 1=Goto, 2=Init, 3=Get Info,
        4=Open CFW Device, 5=Close CFW Device

    ulong cfwParam1 – command specific
    ulong cfwParam2 – “ “
    ushort outLength – “ “
    uchar *outPtr – “ “
    ushort inLength – “ “
    uchar *inPtr – “ “
```

```
}
```

*Results Struct:*

```
typedef struct CFWResults
{
    ushort cfwModel      – See cfwModel above
    ushort cfwPosition   – (see the CFW_POSITION enum)
                        0=Unknown, 1 thru 10 = Position 1 thru 10
    ushort cfwStatus     – (see the CFW_STATUS enum)
                        0=Unknown, 1=Idle, 2= Busy
    ushort cfwError      – (see the CFW_ERROR enum)
                        0=No Error, 1=CFW Busy, 2=Bad Command,
                        3=Calibration Error, 4 = Motor Timeout, 5=Bad Model
    ulong cfwResult1     – command specific
    ulong cfwResult2     – “    “
}
```

### ***CFW Command CFWC\_QUERY***

- Use this command to monitor the progress of the **Goto** sub-command. This command takes no additional parameters in the CFWParams. You would typically do this several times a second after the issuing the Goto command until it reports CFWS\_IDLE in the **cfwStatus** entry of the CFWResults. Additionally filter wheels that can report their current position (all filter wheels except the CFW-6A or CFW-8) have that position reported in **cfwPosition** entry of the CFWResults.

### ***CFW Command CFWC\_GOTO***

- Use this command to start moving the color filter wheel towards a given position. Set the desired position in the **cfwParam1** entry with entries defined by the CFW\_POSITION enum.

### ***CFW Command CFWC\_INIT***

- Use this command to initialize/self-calibrate the color filter wheel. All SBIG color filter wheels self calibrate on power-up and should not require further initialization. We offer this option for users that experience difficulties with their color filter wheels or when changing between the CFW-2 and CFW-5 in the ST-5C/237. This command takes no additional parameters in the CFWParams struct.

### ***CFW Command CFWC\_GET\_INFO***

- This command supports several sub-commands as determined by the **cfwParam1** entry (see the CFW\_GETINFO\_SELECT enum). Command CFWG\_FIRMWARE\_VERSION returns the version of the CFW firmware in the **cfwResults1** entry of the CFWResults and the number of filter positions the CFW supports in the **cfwResults2** entry, commands CFWG\_DATA\_REGISTERS and CFWG\_CAL\_DATA are for internal SBIG use only and all other commands are undefined.

### ***CFW Commands CFWC\_OPEN\_DEVICE and CFWC\_CLOSE\_DEVICE***

- These commands are used to Open and Close any OS based communications port associated with the CFW and should proceed the first command sent and follow the last command sent to the CFW. While strictly only required for the RS-232 version of the CFW-10 calling these commands is a good idea for future compatibility.

For the RS-232 based CFW-10 set the **cfwParam1** entry to one of the settings CFW\_COM\_PORT enum to indicate which PC COM port is used to control the CFW-10. Again, only the RS232 controlled CFW-10 requires these calls.

*Notes:*

- The CFW Command takes pointers to CFWParams as parameters and CFWResults as results.
- Set the **cfwModel** entry in the CFWParams to the type of filter wheel you want to control. The same value is returned in the **cfwModel** entry of the CFWResults. If you select the CFWSEL\_AUTO option the driver will use the most appropriate model.
- The CFW Command is a single API call that supports multiple sub-commands through the **cfwCommand** entry in the CFWParams. Each of the sub-commands requires certain settings of the CFWParams entries and returns varying results in the CFWResults. Each of these sub-commands is discussed in detail above.
- As with all API calls the **CFW** Command returns an error code. If the error code is CE\_CFW\_ERROR, then in addition the **cfwError** entry in the CFWResults further enumerates the error.
- With the CFW-6A model color filter wheel set the **cfwParam1** item to the desired calibrated pulse width in microseconds for the **Goto** and **Init** sub-commands.
- With the CFW-1602 use the **Goto** command to Position 4 to access a special dark position. After issuing this command you must do a **CFW Init** command before going to any other positions.

### 3.4.11 Motor Focus

The Motor Focus Command is a high-level API for controlling SBIG Motor Focus accessories. It supports the new ST Motor Focus unit and will be expanded as required to support new models in the future.

*Parameters Struct:*

```
typedef struct MFParams
{
    enum    mfModel – (see the MF_MODEL_SELECT enum)
             0=Unknown, 1=Auto Select, 2=ST Motor Focuser
    enum    mfCommand – (see the MF_COMMAND enum)
             0=Query, 1=Goto, 2=Init, 3=Get Info, 4=Abort
    long    mfParam1   - command specific
    long    mfParam2   - “ “
    ushort  outLength  - “ “
    uchar   *outPtr    - “ “
    ushort  inLength   - “ “
    uchar   *inPtr     - “ “
}
```

*Results Struct:*

```
typedef struct MFResults
{
    enum mfModel - (see the mfModel above)
    long mfPosition - (position of the Motor Focus, 0=Center, signed)
    enum mfStatus - (see the MF_STATUS enum)
    enum mfError - (see the MF_ERROR enum)
```



```

    long mfResult1 - command specific
    long mfResult2 - “    “
}

```

### ***Motor Focus Command MFC\_QUERY***

- Use this command to monitor the progress of the **Goto** sub-command. This command takes no additional parameters in the MFParams. You would typically do this several times a second after the issuing the Goto command until it reports MFS\_IDLE in the **mfStatus** entry of the MFResults. Motor Focus accessories report their current position in the **mfPosition** entry of the MFResults struct where the position is a signed long with 0 designating the center of motion or the home position. Also the Temperature in hundredths of a degree-C is reported in the **mfResult1** entry.

### ***Motor Focus Command MFC\_GOTO***

- Use this command to start moving the Motor Focus accessory towards a given position. Set the desired position in the **mfParam1** entry. Again, the position is a signed long with 0 representing the center or home position.

### ***Motor Focus Command MFC\_INIT***

- Use this command to initialize/self-calibrate the Motor Focus accessory. This causes the Motor Focus accessory to find the center or Home position. You *can not* count on SBIG Motor Focus accessories to self calibrate upon power-up and should issue this command upon first establishing a link to the Camera. Additionally you should retain the last position of the Motor Focus accessory in a parameter file and after initializing the Motor Focus accessory, you should return it to its last position. Finally, note that this command takes no additional parameters in the MFParams struct.

### ***Motor Focus Command MFC\_GET\_INFO***

- This command supports several sub-commands as determined by the **mfParam1** entry (see the MF\_GETINFO\_SELECT enum). Command MFG\_FIRMWARE\_VERSION returns the version of the Motor Focus firmware in the **mfResults1** entry of the MFResults and the Maximum Extension (plus or minus) that the Motor Focus supports is in the **mfResults2** entry. The MFG\_DATA\_REGISTERS command is internal SBIG use only and all other commands are undefined.

### ***Motor Focus Command MFC\_ABORT***

- Use this command to abort a move in progress from a previous **Goto** command. Note that this will not abort an Init.

#### ***Notes:***

- The **Motor Focus** Command takes pointers to MFParams as parameters and MFResults as results.
- Set the **mfModel** entry in the MFParams to the type of Motor Focus accessory you want to control. The same value is returned in the **mfModel** entry of the MFResults. If you select the MFSEL\_AUTO option the driver will use the most appropriate model and return the model it found in the **mfModel** entry of the MFResults.
- The **Motor Focus** Command is a single API call that supports multiple sub-commands through the **mfCommand** entry in the MFParams. Each of the sub-commands requires certain settings of

the MFParams entries and returns varying results in the MFResults. Each of these sub-commands is discussed in detail above.

- As with all API calls the **Motor Focus** Command returns an error code. If the error code is CE\_MF\_ERROR, then in addition the **mfError** entry in the MFResults further enumerates the error.

### 3.5 General Purpose Commands

The commands discussed in this section are general-purpose commands that do not fall into one of the groups discussed above. They are used by the application to interrogate the driver and camera.

#### 3.5.1 Establish Link

The Establish Link command is used by the application to establish a communications link with the camera. It should be used before any other commands are issued to the camera (excluding the Get Driver Info command).

*Parameters Struct:*

```
struct EstablishLinkParams
{
    ushort sbigUseOnly – leave set to 0, maintained for historical purposes only
}
```

*Results Struct:*

```
struct EstablishLinkResults
{
    enum cameraType - constant specifying the type of camera as specified by the
                        CAMERA_TYPE enum
}
```

*Notes:*

- The EstablishLinkParams struct was modified in version 4 of the driver and no longer specifies the LPT base address. This data is now supplied to the driver through the Open Device command.
- When establishing a link to an ST-237A the **cameraType** is reported as an original ST237\_CAMERA. This was done for maximum compatibility with existing 3<sup>rd</sup> party software packages. The way you distinguish an ST-237A (16 bit A/D) from the ST-237 (12 bit A/D) is by checking the **gain** item from the Get CCD Info command response. If the gain is less than 1.0 (0x100) you are talking to an ST-237A.

#### 3.5.2 Get CCD Info

The Get CCD Info command is used by the application to determine the model of camera being controlled and its capabilities. For future expandability this command allows you to request several types of information. Currently 6 standard requests are supported but as the driver evolves additional requests will be added.

*Parameters Struct:*

```
struct GetCCDInfoParams
{
```

enum **request** - type of CCD information desired (*see the CCD\_INFO\_REQUEST enum*)

- 0 = standard request for Imaging CCD
- 1 = standard request for Tracking CCD
- 2 = extended request for Camera Info
- 3 = extended request for PixCel255/237 Camera Info
- 4 = secondary extended request for Imaging CCD
- 5 = secondary extended request for Tracking CCD
- 6 = extended request for CCD cameras
- 7,8, etc. - reserved for future expansion

}

*Standard Results Struct :*

*requests 0 and 1 -*

struct GetCCDInfoResults0

{

ushort **firmwareVersion** - version of the firmware in the resident microcontroller in BCD format  
(XX.XX, 0x1234 = 12.34)

enum **cameraType** - constant specifying the type of camera, (see CAMERA\_TYPE enum in  
SBIGUDRV.H)

char **name**[64] - null terminated string containing the name of the camera

ushort **readoutModes** - number of readout modes supported

struct **readoutInfo**[20]

{

ushort **mode** - readout mode to pass to the Readout Line command

ushort **width** - width of image in pixels

ushort **height** - height of image in pixels

ushort **gain** - a four digit BCD number specifying the amplifier gain in e<sup>-</sup>/ADU in the  
XX.XX format.

ulong **pixelWidth** - an eight digit BCD number specifying the pixel width in microns in  
the XXXXXX.XX format.

ulong **pixelHeight** - an eight digit BCD number specifying the pixel height in microns in  
the XXXXXX.XX format.

}

}

*request 2 –*

struct GetCCDInfoResults2

{

ushort **badColumns** - number of bad columns in imaging CCD

ushort **columns**[4] - bad columns

enum **imagingABG** - type of Imaging CCD, 0= No ABG Protection, 1 = ABG Present

char **serialNumber**[10] - null terminated serial number string

}

*request 3 - For the PixCel255/237*

```
struct GetCCDInfoResults3
```

```
{
    enum adSize - 0 = Unknown, 1 = 12 bits, 2 = 16 bits
    enum FilterType - 0 = Unknown, 1 = External, 2 = 2 Position, 3 = 5 Position
}
```

*requests 4 and 5 - For all cameras*

```
struct GetCCDInfoResults4
```

```
{
    ushort capabilitiesBits – Set of bits for additional capabilities:
        b0:    0 = CCD is Full Frame Device, 1 = CCD is Frame Transfer Device,
        b1:    0 = No Electronic Shutter, 1 = Interline Imaging CCD with Electronic Shutter and
                millisecond exposure capability
        b2:    0 = No hardware support for external Remote Guide Head, 1 = Detected hardware
                support for external Remote Guide Head.
        b3:    1 = Supports the special Biorad TDI acquisition mode.
        b4:    1 = AO8 detected.
        b5:    1 = Camera contains an internal frame buffer.
        b6:    1 = Camera requires the StartExposure2 command instead of the older deprecated
                StartExposure command.
        Other: See the CB_XXX_XXX definitions in the sbigurdv.h header file.
    ushort dumpExtra – Number of unbinned rows to dump to transfer image area to
                        storage area
}
```

*requests 6 - For all cameras*

```
struct GetCCDInfoResults6
```

```
{
    unsigned long cameraBits – Set of bits for additional capabilities:
        b0:    0 = STX camera, 1 = STXL camera
        b1:    0 = Mechanical shutter, 1 = No mechanical shutter (only an electronic shutter)
        b2 – b31: reserved for future expansion

    unsigned long ccdBits – Set of bits for additional capabilities:
        b0:    0 = Imaging Mono CCD, 1 = Imaging Color CCD
        b1:    0 = Bayer color matrix, 1 = Truesense color matrix
        b2 – b31: reserved for future expansion

    unsigned long extraBits – Set of bits for additional capabilities:
        b0 – b31: reserved for future expansion
}
```

*Notes:*

- The ST-7/8/etc supports types 0, 1, 2, 4 and 5 requests. The PixCel255/237 supports request types 0, 3, 4 and 5.

- The CCD cameras now support **request 6**. Bit b0 of the **cameraBits** defines the STX/STXL camera, bit b1 if set denotes cameras with an electronic shutter only, so no mechanical shutter is installed in this case.  
Bit b0 of the **ccdBits** defines Mono/Color CCD cameras. If this bit is set, ie. Color CCD chip is installed, then bit b1 defines which color matrix (Bayer or Truesense) is bonded to the CCD. If mono CCD is used, than the bit b1 has no meaning.  
An **extraBits** is reserved for future expansion.
- A zero in the height field of the readoutInfo struct signifies the xN mode vertically.
- Mode 9 with 9x9 binning on some sensors (like the ST-1K) will report the maximum pixel height/width of 99.99 even though the binned pixels are actually larger. If you need the correct pixel size multiply the pixel size from mode 0 by the binning factor.
- Requests 4 and 5 are for use when you bypass the **Start Exposure/End Exposure** commands for reading out the CCD such as when you are doing millisecond type exposures with the Tracking CCD and an AO-7. Normally the **End Exposure** command handles transferring the imaging area and dummy rows of these CCDs to the storage area but when you bypass it you must provide for that transfer by adding the number specified in the *dumpExtra* field to the amount of lines you want to discard with the **Dump Lines** command. If you don't supplement the number of lines to dump you will be digitizing data from the storage area of the CCD and won't see any star images at all.
- Request 5 will indicate whether the camera supports the external Remote Guide Head.
- See the notes for the Establish Link command above.
- There are several sub-models of the ST-L based upon the type of CCD they contain but they all report the same value in the **cameraType** field (STL\_CAMERA) of request 0. The name field includes the sub-model.
- Single Shot Color cameras like the ST-2000XCM will have "Color" in the name item. You can use this to decide whether you want to interpret the color information or not. You can also use request 6 (bit b0 of the ccdBits) defined above to distinguish between the mono/color CCD camera.

### 3.5.3 Get Turbo Status

The current driver does not use this command. It was added in a previous version and never removed. It could be reassigned in the future.

### 3.5.4 Query Command Status

The Query Command Status command is used to monitor the progress of a previously requested command. Typically this will be used to monitor the progress of an exposure, relay closure or CFW-6A move command.

*Parameters Struct:*

```
struct QueryCommandStatusParams
{
    ushort command - command of which the status is desired
}
```

*Results Struct:*

```
struct QueryCommandStatusResults
{
    ushort status - command status
}
```

### 3.5.5 Miscellaneous Control

The Miscellaneous Control command is used to control the Fan, LED, and shutter. The camera powers up with the Fan on, the LED on solid, and the shutter closed. The driver flashes the LED at the low rate while the Imaging CCD is integrating, flashes the LED at the high rate while the Tracking CCD is integrating and sets it on solid during the readout.

*Parameters Struct:*

```
struct MiscellaneousControlParams
{
    LOGICAL fanEnable - set TRUE to turn on the Fan
    enum shutterCommand – (see the SHUTTER_COMMAND enum)
        0=leave shutter alone, 1=open shutter, 2=close shutter, 3=reinitialize shutter, 4=open ST-
        L (STX) external shutter, 5=close ST-L (STX) external shutter
    enum ledState – (see the LED_STATE enum)
        0=LED off, 1=LED on, 2=LED blink at low rate, 3=LED blink at high rate
}
```

The status for this command (from the Query Command Status Command) consists of the following bit fields:

- b7-b0 - Shutter edge - This is the position the edge of the shutter was detected at for the last shutter move. Normal values are 7 thru 9. Any other value including 255 indicates a shutter failure and the shutter should be reinitialized.
- b8 - the Fan is enabled when this bit is 1
- b10b9 - Shutter state, 0=open, 1=closed, 2=opening, 3=closing
- b12b11 - LED state, 0=off, 1=on, 2=blink low, 3=blink high

*Notes:*

- The ST-7/8/etc have a shutter, LED and fan but no filter wheel. To position the CFW-6/8 attached to these cameras use the Pulse Out command.
- The PixCel255 has no shutter, fan control or LED but does have a vane/filter wheel. The settings of the fanEnabled, shutterCommand and ledState parameters are ignored and should be set to 0.
- The PixCel237 has no shutter or LED but does have a fan that can be controlled (turned on and off) and vane/filter wheel. The settings of the shutterCommand and ledState parameters are ignored and should be set to 0. The fanEnable is used to control the fan.
- The ST-7/8/etc will cease communications with the PC while the reinitialize shutter command is in progress. After the driver tells the camera to reinitialize the shutter the driver delays 3 seconds for the process to complete.
- For this command the ST-L's external shutter **DOES NOT** mimic the internal shutter (like it does with the **Start Exposure** command). To control the ST-L's external shutter in the Remote Tracking Head pass to **shutterCommand** the SC\_OPEN\_EXT\_SHUTTER and SC\_CLOSE\_EXT\_SHUTTER values.
- On the STX/STT setting the **fanEnable** field to 0 turns off the fan. Setting **fanEnable** to 1 (or greater than 100) sets the fan to Auto-Speed Control where the fan speed is determined by the STX/STT firmware. Setting **fanEnable** to 2 through 100 sets the fan to manual control at 2 to 100% speed. Note that effective manual speed control is achieved with values between 20 and 100 as the fan doesn't really start spinning at manual speeds below 20.

### 3.5.6 Update Clock

The Update Clock no longer supported and calls to this command will do nothing. On older Windows 95/98/Me based systems run or spawn the SETCLOCK.EXE program to restore the system clock for lost time during the readout of parallel port based cameras. Other Windows versions or systems using USB cameras don't lose time.

### 3.5.7 Read Offset

The Read Offset command is used to measure the CCD's offset. In the SBIG cameras the offset is adjusted at the factory and this command is for testing or informational purposes only.

*Parameters Struct:*

```
struct ReadOffsetParams
{
    enum ccd - the CCD to measure offset (see the CCD_REQUEST enum)
        0 = Imaging CCD
        1 = Tracking CCD
        2 = External Tracking CCD in ST-L, STX only
}
```

*Results Struct:*

```
struct ReadOffsetResults
{
    ushort offset - the CCD's offset
}
```

### 3.5.8 Read Offset 2

The Read Offset 2 command is used to measure the CCD's offset and the noise in the readout register. In the SBIG cameras the offset is adjusted at the factory and this command is for testing or informational purposes only.

*Parameters Struct:*

```
struct ReadOffsetParams
{
    enum ccd - the CCD to measure offset (see the CCD_REQUEST enum)
        0 = Imaging CCD
        1 = Tracking CCD
        2 = External Tracking CCD in ST-L, STX only
}
```

*Results Struct:*

```
struct ReadOffsetResults2
{
    ushort offset - the CCD's offset
    double rms - noise in the ccd readout register in ADUs rms
}
```

### 3.5.9 Get US Timer

This command is of extremely limited (and unknown) use. When you have established a link to a parallel port based camera under Windows NT/2000/XP this command returns a counter with 1 microsecond resolution. Under all other circumstances the counter is zero.

### 3.5.10 Set/Get IRQ

This command allows you to control the IRQ priority of the driver under Windows NT/2000/XP. The default settings should work fine for all users and these commands should not need to be used.

We use three settings in our CCDOPS software: High = 27, Medium = 15, Low = 2. Under fast machines Low will work fine. On slower machines the mouse may get sluggish unless you select the Medium or High priority.

### 3.5.11 Get Link Status

This command returns the status of the communications link established with the camera.

*Results Struct:*

```
struct
{
    LOGICAL linkEstablished – TRUE when a link has been established
    ushort baseAddress – base address of the LPT port
    ushort cameraType – CAMERA_TYPE enum
    ulong comTotal – total number of communications with camera
    ulong comFailed – total number of failed communications with camera
}
```

### 3.5.12 Get Error String

This command returns a null terminated C string in English (not Unicode) corresponding to the passed error number. It's handy for reporting driver level errors to the user.

### 3.5.13 Set Driver Control

This command is used to modify the behavior of the driver by changing the settings of one of the driver control parameters. Driver options can be enabled or disabled with this command.

*Parameters Struct:*

```
struct SetDriverControlParams
{
    enum controlParameter - the parameter to modify
                        (see the DRIVER_CONTROL_PARAM enum)
    long controlValue – the value of the control parameter
}
```

*Notes:*

- The DCP\_USB\_FIFO\_ENABLE parameter defaults to TRUE and can be set FALSE to disable the FIFO and associated pipelining in the USB cameras. You would do this for example in applications using Time Delay Integration (TDI) where you don't want data in the CCD digitized until the actual call to ReadoutLine is made.



- The DCP\_CALL\_JOURNAL\_ENABLE parameter defaults to FALSE and can be set to TRUE to have the driver broadcast Driver API calls. These broadcasts are handy as a debug tool for monitoring the sequence of API calls made to the driver. The broadcasts can be received and displayed with the Windows based SBIGUDRVJournalRx.exe application.  
Only use this for testing purposes and do not enable this feature in your released version of your application as the journaling mechanism can introduce minor artifacts in the readout.
- The DCP\_IVTOH\_RATIO parameter sets the number of Vertical Rows that are dumped (fast) before the Horizontal Register is dumped (not as fast) in the DumpRows command for Parallel Port based cameras. This is a very specialized parameter and you should think hard about changing it if you do. The default of 5 for the IHTOV\_RATIO has been determined to offer a good compromise between the time it takes to clear the CCD or Dump Rows and the ability to effectively clear the CCD after imaging a bright object. Finally should you find it necessary to change it read the current setting and restore it when you're done.
- The DCP\_USB\_FIFO\_SIZE parameter sets the size of the FIFO used to receive data from USB cameras. The default and maximum value of 16384 yields the highest download speeds. Lowering the value will cause the camera to digitize and download pixels in smaller chunks. Again this is a specialized parameter that 99.9% of programs out there will have no need for changing.
- The DCP\_USB\_PIXEL\_DL\_ENABLE parameter allows disabling the actual downloading of pixel data from the camera for testing purposes. This parameter defaults to TRUE.
- The DCP\_HIGH\_THROUGHPUT parameter allows configuring the driver for the highest possible imaging throughput at the expense of image noise and or artifacts. This parameter defaults to FALSE and you should only enable this for short periods of time. You might use this in Focus mode for example to get higher image throughput but you should never use it when you are taking keeper images. It does things that avoid timed delays in the camera like leaving the shutter motor on all the time, etc. At this time this feature is supported in the driver but not all cameras show a benefit from its use.
- The DCP\_VDD\_OPTIMIZED parameter defaults to TRUE which lowers the CCD's Vdd (which reduces amplifier glow) only for images 3 seconds and longer. This was done to increase the image throughput for short exposures as raising and lowering Vdd takes 100s of milliseconds. The lowering and subsequent raising of Vdd delays the image readout slightly which causes short exposures to have a different bias structure than long exposures. Setting this parameter to FALSE stops the short exposure optimization from occurring.
- The DCP\_AUTO\_AD\_GAIN parameter defaults to TRUE whereby the driver is responsible for setting the A/D gain in USB cameras. Setting this to FALSE allows overriding the driver imposed A/D gains.
- The DCP\_NO\_HCLKS\_FOR\_INTEGRATION parameter defaults to FALSE and setting it to TRUE disables the horizontal clocks during exposure integration and is intended for SBIG testing only.
- The DCP\_TDI\_MODE\_ENABLE parameter defaults to FALSE and setting it to TRUE enables the special Biorad TDI mode.
- The DCP\_VERT\_FLUSH\_CONTROL\_ENABLE parameter defaults to TRUE and setting it to FALSE it disables the background flushing of the vertical clocks of KAI CCDs during exposure integration and is intended for SBIG testing only.
- The DCP\_ETHERNET\_PIPELINE\_ENABLE parameter defaults to FALSE and setting it to TRUE can increase the throughput of Ethernet based cameras like the STX & STT but doing so is not recommended for robust operation.

- The DCP\_FAST\_LINK parameter defaults to FALSE and setting it to TRUE speeds up the Establish Link command by not dumping the pixel FIFOs in the camera. It is used internally to speed up the Query USB and Query Ethernet commands.
- The DCP\_COLUMN\_REPAIR\_ENABLE defaults to FALSE and setting it to TRUE causes the Universal Driver Library to repair up to 7 columns in the Imaging CCD automatically. This is done in conjunction with column data stored in nonvolatile memory in the cameras. Under Windows the setting of this parameter persists in the Registry through the setting of the *HKEY\_CURRENT\_USER\Software\SBIG\SBIGUDRV\Filter\ColumnRepairEnable* setting.
- The DCP\_WARM\_PIXEL\_REPAIR\_ENABLE defaults to Zero and setting it to 1 through 8 causes the Universal Driver Library to repair warm pixels in the Imaging CCD automatically. A setting of 8 replaces approximately 5% of pixels and a setting of 1 replaces approximately 1 in a million. A decrease of 1 in the setting replaces approximately 1/10<sup>th</sup> the number of pixels of the higher setting (7  $\approx$  0.5%, 6  $\approx$  0.05%, etc). Under Windows the setting of this parameter persists in the Registry through the setting of the *HKEY\_CURRENT\_USER\Software\SBIG\SBIGUDRV\Filter\WarmPixelRepairEnable* setting.
- The DCP\_WARM\_PIXEL\_REPAIR\_COUNT parameter returns the total number of pixels replaced in the last image by the Warm Pixel Repair routine described above. You can use this parameter to tweak the DCP\_WARM\_PIXEL\_REPAIR\_ENABLE parameter to filter as many warm pixels as your application requires.

*Notes:*

- There is one set of parameters for the whole DLL vs. one per handle.

### 3.5.14 Get Driver Control

This command is used to query the setting of one of the driver control parameters.

*Parameters Struct:*

```
struct GetDriverControlParams
{
    enum controlParameter - the parameter to modify
                        (see the DRIVER_CONTROL_PARAM enum)
}
```

*Results Struct:*

```
struct GetDriverControlResults
{
    long controlValue – the value of the control parameter
}
```

*Notes:*

- See the Set Driver Control command above.
- The DCP\_WARM\_PIXEL\_REPAIR\_COUNT parameter returns the number of pixels replaced in the previous image with the STF-8300's Auto Filter capability.

### 3.5.15 USB AD Control

This command is used to modify the USB cameras A/D gain and offset registers.

*Parameters Struct:*

```

struct USBADControlParams
{
    enum command – Imaging or Tracking CCD Gain or Offset
                    (see the USB_AD_CONTROL_COMMAND enum)
    short data – command specific
}

```

*Notes:*

- This command is intended for OEM use only. The typical application does not need to use this command as the USB cameras initialize the A/D to factory set defaults when the camera powers up.
- For the *USB\_AD\_IMAGING\_GAIN* and *AD\_USB\_TRACKING\_GAIN* **commands** the allowed setting for the **data** parameter is 0 through 63. The actual Gain of the A/D (in Volts/Volt) ranges from 1.0 to 6.0 and is determined by the following formula:

$$\text{Gain} = 6.0 / ( 1.0 + 5.0 * ( (63 - \text{data}) / 63 ) )$$

Note that the default A/D Gain set by the camera at power up is 1.2 for the Imaging CCD and 2.0 for the Tracking CCD. Furthermore, the gain item reported by the **Get CCD Info** command will always report the default factory-set gain and will not change based upon changes made to the A/D gain by this command.

- For the *USB\_AD\_IMAGING\_OFFSET* and *USB\_AD\_TRACKING\_OFFSET* **commands** the allowed setting for the **data** parameter is –255 through 255. Positive offsets increase the video black level in ADUs. The cameras are programmed at the factory to typically have a 900 to 1000 ADU black level offset.

**3.5.16 Query USB**

This command is used to query the USB bus and detect up to eight cameras. This allows the user to have multiple USB cameras and for the programmer to connect with specific cameras. This command takes no parameters and returns data for up to eight cameras.

*Results Struct:*

```

typedef struct QueryUSBResults
{
    unsigned short camerasFound – total number of cameras found, 0 through 8
    struct QUERY_USB_INFO usbInfo[8] – data for up to 8 cameras
        MY_LOGICAL cameraFound – TRUE if camera found
        enum CAMERA_TYPE cameraType – type of camera found
        char name[64] – name of camera (like from GetCCDInfo command)
        char serialNumber[10] – camera's serial number (like from GetCCDInfo command)
}

```

*Notes:*

- To Establish a link to a specific camera specify DEV\_USB1, DEV\_USB2 ... DEV\_USB8 in the device field of the **Open Device** command. DEV\_USB1 corresponds to the camera described in usbInfo[0], DEV\_USB2 corresponds to usbInfo[1], etc. If you specify USB\_DEV in Open Device it opens the next available device.

- You should call this command after calling **Open Driver** but before calling **Open Device**.

### 3.5.17 Query Ethernet

This command is used to query and detect up to eight cameras Ethernet based cameras. This allows the user to have multiple Ethernet cameras and for the programmer to connect with specific cameras. This command takes no parameters and returns data for up to eight cameras.

*Results Struct:*

```
typedef struct QueryEthernetResults
{
    unsigned short camerasFound – total number of cameras found, 0 through 8
    struct QUERY_ETHERNET_INFO ethernetInfo[8] – data for up to 8 cameras
        MY_LOGICAL cameraFound – TRUE if camera found
        ulong ipAddress – IP address of camera
        enum CAMERA_TYPE cameraType – type of camera found
        char name[64] – name of camera (like from GetCCDInfo command)
        char serialNumber[10] – camera's serial number (like from GetCCDInfo command)
}
```

*Notes:*

- To establish a link to a specific camera specify the discovered IP Address in the IP Address field of the **Open Device** command.
- You should call this command after calling **Open Driver** but before calling **Open Device**.

### 3.5.18 Get Pentium Cycle Count

This command is used to read a Pentium processor's internal cycle counter. Pentium processors have a 32 or 64 bit register that increments every clock cycle. For example on a 1 GHz Pentium the counter advances 1 billion counts per second. This command can be used to retrieve that counter.

*Parameters Struct:*

```
typedef struct GetPentiumCycleCountParams
{
    short rightShift – number of bits to shift the results to the right (dividing by 2)
}
```

*Results Struct:*

```
typedef struct GetPentiumCycleCountResults
{
    unsigned long countLow – lower 32 bits of the Pentium cycle counter
    unsigned long countHigh – upper 32 bits of the Pentium cycle counter
}
```

*Notes:*

- Only call this function if you are sure your code is running on a Pentium processor.
- As the Pentium counter increments every clock cycle and the rate depends on the CPU clock rate, if you want to use this command for generating precise delays you'll need to calibrate the counter over a period of time like 1 or 10 seconds.

- As the Pentium clock rates are quite high (typically hundreds of Megahertz or higher) this counter has much higher resolution than most programs need. Furthermore the 64-bit result can be cumbersome. Setting the **rightShift** item in the parameters allows reducing the resolution. For example on a 1 GHz Pentium setting the rightShift to 10 divides the results by  $2^{10}$  or 1024 returning a good approximation to a 1 microsecond timer that overflows the **countLow** only once every 4295 seconds.

### 3.5.19 RW USB I2C

This command is used read or write data to the USB cameras I2C expansion port.

*Parameters Struct:*

```
typedef struct RWUSBI2CParams
{
    unsigned char address – Address to read from or write to
    unsigned char data – Data to write to the external I2C device, ignored for read
    MY_LOGICAL write – TRUE when write is desired , FALSE when read is desired
    unsigned char deviceAddress – Device Address of the I2C peripheral
}
```

*Results Struct:*

```
typedef struct RWUSBI2CResults
{
    unsigned char data – Data read from the external I2C device
}
```

*Notes:*

- This command is typically called by SBIG code in the Universal Driver. If you think you have some reason to call this function you should check with SBIG first.

### 3.5.20 Bit IO

This command is used read or write control bits in the USB cameras.

*Parameters Struct:*

```
typedef struct BitIOParams
{
    unsigned short bitOperation – 0=Write, 1=Read
    unsigned short bitName – 0=Read Power Supply Low Voltage, 1=Write Genl. Purp. Bit 1,
                           2=Write Genl. Purp. Bit 2, 3=Read Genl. Purp. Bit 3
    MY_LOGICAL setBit - 1=Set Bit, 0=Clear Bit
}
```

*Results Struct:*

```
typedef struct BitIOResults
{
    MY_LOGICAL bitIsSet - 1=Bit is set, 0=Bit is clear
}
```

*Notes:*

- On the ST-L camera you can use this command to monitor whether the input power supply has dropped to the point where you ought to warn the user. Do this by issuing a Read operation on bit 0 and if that bit is set the power has dropped below 10 Volts.

### 3.5.21 Column EEPROM

This command is used read or write the STF-8300's Column Repair data stored in the camera for use with that camera's Auto Filter Capability.

*Parameters and Results Struct:*

```
typedef struct
{
    MY_LOGICAL writeData – Set to 0 to read the data from the camera and 1 to write it
    unsigned short columns[7] – Specify up to 7 columns to repair
    unsigned short flags – not used at this time
}
ColumnEEPROMParams, ColumnEEPROMResults;
```

*Notes:*

- The left most column is column 1 (not zero). Specifying a column zero doesn't filter any columns.
- This command is somewhat unique in that the *Parameters* and the *Results* are the same struct.
- To enable column filtering you must use this command and also the Set Driver Control command to set the DCP\_COLUMN\_REPAIR parameter to 1.

### 3.5.22 Customer Options

This command is used read or write the STX/STXL/STT's customer options which are:

- images includes overscan regions
- enable window heater
- preflash CCD
- Vdd normally off

*Parameters and Results Struct:*

```
typedef struct
{
    MY_LOGICAL bSetCustomerOptions - set to 0 to read the data and 1 to write it
    MY_LOGICAL bOverscanRegions
    MY_LOGICAL bWindowHeater
    MY_LOGICAL bPreflashCcd
    MY_LOGICAL bVddNormallyOff
}
CustomerOptionsParams, CustomerOptionsResults;
```

*Notes:*

- The `bSetCustomerOptions` of the `CustomerOptionsParams` structure is always copied to the `CustomerOptionsResults` structure on return.
- To read customer options from the camera, just set the `bSetCustomerOptions` to zero inside the `CustomerOptionsParams` structure. Other fields in this structure are not used. The `CustomerOptionsResults` returns requested values.
- To write customer options to the camera, set `bSetCustomerOptions` of the `CustomerOptionsParams` to one, including all requested settings. Value 1 means you want to set corresponding value inside the camera, value 0 means you want to clear that value. The `CustomerOptionsResults` returns copies of all input parameters.

### 3.5.23 Debug Log

This command is used to enable / disable writing messages to a log file.

The type of log messages can be enabled or disabled using corresponding bits as follows:

Debug Log Bit Flags monitors:

- `DLF_DRV` = 0x0001 - basic driver commands
- `DLF_CFW` = 0x0002 - CFW commands
- `DLF_AO` = 0x0004 - AO commands
- `DLF_READOUT_LINE` = 0x0008 - CCD Readout commands
- `DLF_TEMPERATURE` = 0x0010 - Temperature commands
- `DLF_ETHERNET` = 0x0020 - Ethernet commands
- `DLF_USB` = 0x0040 - USB commands
- `DLF_MICRO_COMMAND` = 0x0080 - driver's Micro commands
- `DLF_READ_COM_PIPE` = 0x0100 - Read Com pipe commands
- `DLF_WRITE_COM_PIPE` = 0x0200 - Write Com pipe commands
- `DLF_READ_PIXEL_PIPE` = 0x0400 - Read Pixel pipe commands
- `DLF_READ_ALT_PIPE` = 0x0800 - Read Alternate pipe commands
- `DLF_WRITE_ALT_PIPE` = 0x1000 - Write Alternate pipe commands
- `DLF_QUERY_CMD_STATUS` = 0x2000 - Query Command Status

*Parameters Struct:*

```
typedef struct
{
    unsigned short debugLogFlags;
}
DebugLogParams;
```

*Notes:*

- There is one set of parameters for the whole DLL vs. one per handle.

- Right now, this function is implemented in Windows / Mac OSX operating systems only.  
The log files are: Windows: C:\sbigudrvlog.txt, OSX: /Library/Logs/sbigudrv.log



## 4. Windows Based Utility Programs

This section describes several Windows based tool programs that SBIG provides to help you get your custom application up and running.

### 4.1 SBIGDriverChecker.exe and SBIGDriverChecker64.exe

These programs check the drivers installed on the system against the “SBIG Driver” directory. The SBIGDriverChecker.exe is for 32 bit versions of Windows and the SBIGDriverChecker64.exe is for 64 bit version of Windows but will run on 32 bit versions of Windows XP and later. It reports whether the drivers installed are current and allows you to download the latest drivers from SBIG’s Servers with the Download button then update all drivers with the Update button. Please read the “Installing USB.pdf” and note that you should distribute this program and the SBIG Drivers directory with your code so users can install the latest drivers. The latest version and the associated drivers can be downloaded from the Software Downloads page of our home page.

This utility needs to be run in Administrator mode and depends on the SBIG Drivers directory being in the same directory as the utility. Also note that if you try to run this program from your installer and you get the error message “Could not find the SBIG Drivers” it’s because the Installer is not setting the current working directory to the utility directory. You can get around this if you pass this program a command line argument with the full path to the utility. For example if these items are installed into the “C:\Program Files\My Company\My Program” directory pass the SBIGDriverChecker.exe utility “C:\Program Files\My Company\My Program\SBIGDriverChecker.exe” without the quotes as the command line argument.

The program usually requires user intervention to update the drivers, which we feel is the safest way to go. You can however pass the program the “/s” command line argument in which case it will run in “silent” mode and automatically update the drivers and then exit.

### 4.2 EthSim2.exe

The EthSim2.exe tool allows you to simulate various model cameras on an Ethernet network. This is handy for testing your programs without a camera. To use EthSim do the following:

- Copy the EthSim2.exe program and your Test Program to computers connected together through an Ethernet LAN. The two programs can run on two different computers or the same computer so long as that computer has Ethernet. You can use CCDOPS as the Test Program.
- Double-click the EthSim2 icon or run EthSim2.exe from the command line. EthSim2 takes an optional command line parameter specifying the type of camera to simulate. For example, to simulate an ST-2K you would run: ethsim2 2K  
with no command line parameter it defaults to an ST-7.
- Select Ethernet for the Comm link and enter the IP address of the computer running EthSim2.exe
- Establish a link to the simulated camera

EthSim2 will generate random pixel data for both the Imaging and Tracking CCDs and simulates the following types of cameras:

ST-5C, ST-237, ST-237A, ST-7, ST-8, ST-9, ST-10, ST-1K, ST-2K, ST-4K, the STL family, ST-402, ST-1602, ST-3200, ST-8300, STX and the STT family.

### 4.3 SBIGUDRVJournalRx.exe

The SBIGUDRVJournalRx.exe tool is used in conjunction with the DCP\_CALL\_JOURNAL\_ENABLE Driver Control Parameter of the **Set Driver Control** command. When that parameter is set TRUE the driver broadcasts API calls. SBIGUDRVJournalRx.exe receives and displays those broadcasts. The broadcasts essentially tell you the sequence of API calls made to the driver.

SBIGUDRVJournalRx has a scrolling window that shows the API call sequence. In addition the program has the following controls:

- **Clear** – Clicking this button clears the log window
- **Filters** – Clicking this button brings up a dialog that allows you to filter out any or all of the API calls from the displayed log window.
- **Show Log** – Uncheck this button to disable the display of API calls
- **Send Msg** – Clicking this button simulates a broadcast from the driver to the program using the WParam and LParam data entry boxes. Each time you click this button a message should appear in the log.

Finally note that if you are not receiving API broadcasts messages from the driver it is probably because you did not enable journaling with the Set Driver Control command. Just for reference CCDOPS does not enable them.

### 4.4 SetClock.exe

This 16-Bit windows program, which you can distribute with your Application if you find it useful, reads the CMOS clock on the PC and sets the Windows System time to that value. It can take up to 1 second to synchronize the clocks as it waits till the CMOS clock hits the second mark. Typically you would spawn this program at the end of the readout on a Parallel camera on Windows 95/98/Me based systems to adjust for the seconds lost while interrupts were disabled for the readout.

### 4.5 GetPortD.exe

This 16-Bit DOS program is used to determine the I/O address of the LPT ports on Windows 95/98/Me based systems. It reads the LPT address information from DOS, displays the address of the three LPT ports and creates a binary file named PORTADDR.DAT containing those addresses. We use this to allow the user to select “LPT1, LPT2 or LPT3” instead of the I/O address for the Parallel port. Spawn this program, then open and read three unsigned shorts from the PORTADDR.DAT file. The first is the I/O address of LPT1, the second LPT2 and the third LPT3. You then pass this information to the Open Device command.

## 5. Supporting New Cameras and Accessories

This section gives specific instructions about how to support relatively new model cameras and accessories starting with the ST-L.

### 5.1 Supporting the ST-L

The ST-L (Large Format Camera) is a USB 1.1 based camera much like the USB version of the ST-7/8/... with the following significant differences:

- There are several models of the ST-L with various imaging CCDs. These include the 11 Megapixel ST-L-11K with the Kodak KAI-11000 and various other models. While there is only a single **cameraType** (STL\_CAMERA) returned by the **Establish Link** command, the **Get CCD Info** command returns the appropriate pixel information (pixel size and number of pixels) and the **name** includes the model type.
- The camera has an internal 5-position filter wheel referred to as the CFW-L. The carousel allows using 2-inch diameter filter substrates or filters mounted in 48mm cells. To support the CFW-L use the new **CFW** Command described in Section 3.4.10. Unlike the CFW-8, the CFW-L actually detects erroneous moves and reports those back via the **CFW** command. In addition the CFW-L is self-calibrating on power-up but can be recalibrated at any time with the **CFW** command.
- The ST-L firmware does not support generating pulses out the CFW port. The Universal Driver will redirect **Pulse Out** commands to the internal CFW-L making the CFW-L mimic a CFW-8. This allows CFW-8 software to work with the CFW-L but you really should add support for the CFW-L with the **CFW** command described in Section 3.4.10.
- Like the USB ST-7/8/... the ST-L contains an internal tracking CCD. At this time the internal tracking CCD is a TC-237 but you should depend on the results of the **Get CCD Info** command when deciding how to handle the tracking CCD.
- The ST-L also supports an optional plug-in External Guider. The external guider is a cooled TC-237 based camera with shutter for taking dark frames. The pixel details of the external guider are the same as the internal tracking CCD and are returned by the **Get CCD Info** command. To utilize the external guider use the CCD\_EXT\_TRACKING enum value in the **ccd** entry of the command like **Start Exposure**, **End Exposure**, **Readout Line**, etc.
- Since the External Guider is no longer in a fixed orientation relative to the imaging CCD you will need to maintain two separate sets of Tracking Calibration Vectors if your software supports guiding.
- The External Guider's cooling is open loop and is either on or off. It mimics the main cooling in that it is off only when the main cooler is at off. Otherwise it is on. You don't need any special command for dealing with the external cooling.
- The External Guider's shutter mimics the internal shutter with respect to the **Start Exposure** command. When the internal shutter is asked to open the external shutter opens, etc. If however you are using the **Miscellaneous Control** command for shutter control there are separate commands for the internal and external shutters.
- The ST-L runs off a single external 12V power supply. The ST-L's internal power supply monitors the voltage level into the camera and has several stages of warning and protection as that voltage gets reduced. Below 11 Volts a warning light on the side of the camera comes on. Below 10 volts a second warning light comes on. This condition can be detected in software using the **Bit IO** command described in Section 3.5.20. Below 9 volts a third warning light comes on and power to the TE cooler is limited. You should monitor the voltage periodically

with the **Bit IO** command and warn the user when there is a transition from above 10V to below 10V.

## 5.2 Supporting the Single Shot Color Cameras

SBIG has recently introduced the model ST-2001XCM camera, which unlike all our previous cameras, contains a color CCD. With Color CCDs each group of four pixels contains 2 Green pixels, a Red pixel and a Blue pixel arranged in a matrix as show below:



In the literature this is called a “Bayer Pattern” or “Color Filter Array” (CFA). You can tell whether a particular camera contains a color CCD by looking for “Color” in the **name** item returned by the **Get CCD Info** command.

As you readout the data from a color CCD using the **Readout Line** command, even rows (0, 2, 4) start with a Blue pixel and alternate between Blue and Green. Odd rows (1, 3, 5) start with a Green pixel and alternate between Green and Red. For sub-frame readout, if you start at an even row and even column this will still be true. Otherwise the matrix will be offset by one horizontally or vertically.

To make a full resolution RGB image you can interpolate pixels. For example to calculate a Green value at the position of a Blue or Red pixel you could average the four Green pixels above, below, to the right and to the left. Calculating a Red or Blue value at a Green position could be done using an average of the Left and Right pixels, etc. Simple algorithms like this seem to work well for Astrophotos but tend to show a “zipper effect” along color edges. If you want to read about other color extraction algorithms search the web for “Bayer Pattern” or “CFA Demosaicing”.

Finally extracting a monochrome or Luminance image is easy if you apply to the image the 3x3 kernel filter shown below:

1	2	1
2	4	2
1	2	1

16

This yields the equation:  $L = (R + 2G + B) / 4$ .

## 5.3 Supporting the ST-402 Camera and CFW-402 Filter Wheel

In September of 2004 SBIG introduced the ST-402 camera to replace the out of production ST-237A and ST-5C camera models. The ST-402 is a small (4 x 5 x 2.5 inches) single-CCD camera based on the KAF-0402ME CCD with 765 x 510 9-micron square pixels. Like the ST-237/ST-5C the camera has an optional internal 4-position Color Filter Wheel (CFW-402) with Red, Green, Blue and Clear filters. Additionally, the camera has a USB 2 interface (backwards compatible with USB 1.1) that offers downloads rates 3 times as fast as the USB 1.1 versions of the ST-7.

Supporting the ST-402 in software is quite simple and works like every other USB camera in that regard. It acts just like an ST-7 with only an Imaging CCD and with an optional CFW-402 Color Filter Wheel. As with the ST-7, shutter control is through the **Start Exposure** command. Additionally the CFW-402 is supported by the **CFW** command with the **cfwModel** set to **CFWSEL\_CFW402**. Additionally the Red,

Green, Blue and Clear filters being positions are represented by **CFWP\_1** through **CFWP\_4** respectively.

#### 5.4 Supporting the CFW-10 Color Filter Wheel

In September of 2004 SBIG introduced the CFW-10 Color Filter Wheel. The CFW-10 holds up to ten 1 ¼ inch filters in threaded Filter Cells. The CFW-10 interfaces to the USB ST-7 series cameras through those camera's I2C-Aux port, drawing power and control from that port. Optionally the CFW-10 can be used with older Parallel port based cameras or non-SBIG cameras with an external 12V power supply and control via RS-232.

When attached to the USB ST-7 series of cameras the CFW-10 can be controlled with the Universal Driver's **CFW** command the **cfwModel** set to **CFWSEL\_CFW10**. For RS-232 based control of the CFW-10 set the **cfwModel** to **CFWSEL\_CFW10\_SERIAL** and don't forget to call the **CFWC\_OPEN\_DEVICE** and **CFW\_CLOSE\_DEVICE** commands to open and close the PC's COM port. Finally, filter positions 1 through 10 are represented by the **CFWP\_1** through **CFWP\_10** constants.

#### 5.5 Supporting the I<sup>2</sup>C based AO Accessories

Supporting the I<sup>2</sup>C based AO-L and Low-Cost AO-8 has been made transparent to AO-7 users by having the driver automatically detect the type of AO attached to the camera upon establishing a link and then routing the **AO Tip Tilt** commands to the correct port. The only thing that's different from a programmer's standpoint is that the I<sup>2</sup>C AO's can be re-centered at any time with the new **AO Center** command. For further information on the AO-8 see section 5.8 below.

#### 5.6 Supporting the I<sup>2</sup>C based CFW-9, CFW-L8 and CFW-L8G Color Filter Wheels

Supporting the new I<sup>2</sup>C based color filter wheels is simply a matter of using the **CFW** command like every with other filter wheel. You can manually or auto-select the CFW model in the **CFW** command

#### 5.7 Supporting the I<sup>2</sup>C based Motor Focus Accessory

Supporting the new I<sup>2</sup>C based motor focus accessory is accomplished through the new **Motor Focus** command. You can manually or auto-select the Motor Focus model in the **Motor Focus** command

#### 5.8 Supporting the I<sup>2</sup>C based low-cost AO-8

In September of 2007 we introduced the AO-8 which is an I<sup>2</sup>C based replacement for AO-7 for use with the small format ST series of cameras. The AO-8 displaces the image by tilting a transmissive glass plate where as the AO-7 displaced the image by tilting a 90° reflective mirror. This allowed us to make the AO-8 thinner. Here's what you need to know about the AO-8 compared to the older AO-7:

- The SBIG Universal Driver Library detects whether an AO-8 is attached to the camera in the **CC\_ESTABLISH\_LINK** call, and if so it formats and sends AO commands to the I<sup>2</sup>C port. If no AO-8 is detected then AO commands are formatted for an AO-7 and sent to the AO/CFW/SCOPE port. This is completely transparent to any Application level software and only requires that the AO-8 be attached at the time you establish the communications link.
- Compared to an AO-7, the AO-8 calibration vectors are very perpendicular and in fact, due to the transmissive design the Calibration Vectors are actually fixed: *they repeat from unit to unit and camera to camera*. Because of this users may skip the Calibration step if you proved them with default Calibration Vectors as described below:

Relay ->	X+	X-	Y+	Y-
Calibration Vector Standard Orientation	Right	Left	Down	Up
Calibration Vector Rotated Orientation	Left	Right	Up	Down
Speed with 237 Guider	2.25	Same	Same	Same
Speed with 211 Guider	1.0	Same	1.2	Same

Note that the Calibration Vectors change orientation depending on whether the AO-8 is attached in the Standard Orientation with its cable pointed towards the camera's Power Cable or rotated 180° as is required when used with a CFW.

- The **capabilitiesBits** field of the CCD\_INFO\_EXTENDED2\_IMAGING request of the CC\_GET\_CCD\_INFO command contains a field that allows you to find out whether an AO-8 was detected when the link was established to the camera. You can use this to preload the calibration vectors for the AO-8 as described in paragraph 2 above.

### 5.9 Supporting the ST-4000XCM Single Shot Color Camera

In September of 2007 we introduced the ST-4000XCM. The ST-4000XCM is a 4-megapixel single-shot-color camera in the smaller ST-series camera body. Here's what you need to know that's different about the ST-4000XCM:

- When you establish a link to the ST-4000XCM the driver returns a **cameraType** code of ST4K\_CAMERA (18) as defined in the sbigudrv.h file.
- Like the STL version of this camera, the ST-4000XCM has the imager rotated 90° clockwise relative to the tracking CCD. This tends to easily disorient users who are used to the standard orientation. For this reason we recommend *you automatically rotate ST-4000XCM images 90° clockwise* for the user. That's what we do in CCDOps. For full-frame images that's fairly simple because the imager itself is square at 2048 x 2048 pixels. It gets a little more complicated with rectangular sub-frames but it's not impossible. Just remember that displayed images have been rotated 90° CW so rotate any user selected sub-frame coordinates 90° CCW to get to the CCD coordinates.. Usually we like to take care of issues like this in the driver so that they are taken care for you automatically but in this case we could not because the API to the driver is row based, not area based.

### 5.10 Supporting the STX Cameras

Towards the end of 2009 we introduced the STX family of cameras. This line of cameras supports multiple models of CCDs and offers both USB and Ethernet connectivity. While we tried to minimize the impact on developers of adding support for this camera to existing software, the architectural change of including an on-camera frame buffer requires adding the following capabilities to existing software:

- With the STX you need to use the **Start Exposure 2** command instead of the **Start Exposure** command so you can specify the pixel coordinates of the desired readout area. You can use the **Start Exposure 2** command with the older cameras as well.
- Even though you use the **Start Exposure 2** command you should continue to use the **Start Readout** and **End Readout** commands and pass the same values for the binning and image coordinates to these other commands. You should also use the **Dump Lines** command to discard unused rows at the top of the image.

- Please check out the modifications to the **Query Temperature Status** command and the new **Query Ethernet** command. These offer ways to improve support for SBIG cameras and in particular the STX.

### 5.11 Supporting the ST-i Cameras

Towards the end of 2010 we introduced the ST-i camera for shipment in the 2<sup>nd</sup> Quarter of 2011. The ST-i is a low cost, eyepiece sized, uncooled, single CCD camera based on the KAI-0340 CCD with approximately 640 x 480 pixels. The camera has an internal shutter and a standard SBIG tracking port and is available in both monochrome and single shot color versions. There are a couple of special requirements for the ST-i:

- You must use the **Start Exposure 2** command instead of the now deprecated **Start Exposure** command.
- You detect the single shot color version by the presence of the word “Color” in the **name** field as returned by the **Get CCD Info** command just like other color cameras.
- Also as indicated by the **Get CCD Info** command the ST-I supports millisecond exposures.
- The ST-i only supports four readout modes: High Res (1x1), Medium Res (2x2) and two By-N modes (1xN, 2xN). There is no Low Res mode (3x3) like on many of the other cameras. This information is provided by the **Get CCD Info** command.

### 5.12 Supporting the STF-8300 Cameras

In the fall of 2011 we introduced the STF-8300 camera which is a faster readout version of the ST-8300 Camera with the addition of a Full Frame Image Buffer. The changes required to support the STF-8300 include:

- You must use the **Start Exposure 2** command instead of the now deprecated **Start Exposure** command.
- The STF-8300 supports an Auto Filter capability whereby you can configure the camera and Universal Driver Library to automatically repair bad columns and repair warm pixels in the images acquired. See sections 3.5.13, 3.5.14 and 3.5.21 for how to configure the Auto Filter capability.
- The STF-8300 has a fast readout mode where it allows simultaneous pixel digitization and download at the expense of increased low-level artifacts in the images. This effectively cuts the image download time in half. See the description of the Start Exposure 2 command in section 3.2.1.
- The STF-8300 camera supports an new on-chip NxN binning mode where N is 1 to 256. See the Start Exposure 2 command for more details.

### 5.13 Supporting the STT Cameras

In April 2012 we introduced the STT family of cameras. This line of cameras supports multiple models of CCDs and offers both USB and Ethernet connectivity. This camera is similar to the STX family with the exception that only has an internal tracking CCD. There is no support for the Remote Tracking Head. While we tried to minimize the impact on developers of adding support for this camera to existing software, the architectural change of including an on-camera frame buffer requires adding the following capabilities to existing software:

- With the STT you need to use the **Start Exposure 2** command instead of the **Start Exposure** command so you can specify the pixel coordinates of the desired readout area. You can use the **Start Exposure 2** command with the older cameras as well.
- Even though you use the **Start Exposure 2** command you should continue to use the **Start Readout** and **End Readout** commands and pass the same values for the binning and image coordinates to these other commands. You should also use the **Dump Lines** command to discard unused rows at the top of the image.
- Please check out the modifications to the **Query Temperature Status** command and the new **Query Ethernet** command. These offer ways to improve support for SBIG cameras and in particular the STT.

## 6. Revision History

This section details the recent changes to this specification and supporting software since the initial release of the ST-7 Driver.

### Changes Incorporated in Version 1.90

- Added support for the PixCel255.
- Added a response to the EstablishLink command that reports the type of camera found.
- Renamed the main function SBIGUDrvCommand() from ST7Command() since it now supports the PixCel255 in addition to the ST-7/8.
- Added the new Error Code Unknown Camera.
- Added several ABG rates to the abgState item in the Start Exposure command.
- Added status results to the Pulse Out command that reports the current filter in position in the PixCel255.
- Added a type 3 request to the Get CCD Info command to report data on the PixCel255 vane/filter wheel configuration.

### Changes Incorporated in Version 1.96

- Renamed ST-5C to PixCel255.
- Set PixCel255 e-/ADU based upon production hardware.

### Changes Incorporated in Version 2.1

- Made 32 bit version available for Windows 95.
- Added the Open Driver and Close Driver commands.
- Added Driver Not Found, Driver Not Open and Driver Not Closed error codes for Open Driver and Close Driver functions.
- Added the Nx1, Nx2 and Nx3 readout modes for the ST-7/8
- Added the 1x1, 2x2 and 3x3 off-chip binning modes for the ST-7/8

### Changes incorporated through Version 2.6

- Added support for the PixCel237 and the AO-7.

### Changes incorporated through Version 2.70

- Added the **AO Delay** command for generating millisecond level delays.
- Added the **End Readout** command for preparing the CCD for the idle state.
- Added the ability to Freeze the TE Cooler for readout through the **Set Temperature Regulation** command.
- Added the REGULATION\_FOZEN\_MASK bit to the *enabled* field of the **Query Temperature Status** command results for detecting whether the TE cooler is frozen.



## Changes incorporated through Version 3.3

- Added the **Open Device** and **Close Device** commands.
- Added support for the ST-9
- Added support for the ST-10
- Added support for the ST-1K
- Added support for an OEM version of the STV with Parallel Port
- Added support for the ST-237A

## Changes incorporated through Version 4.0

- Added the START\_SKIP\_VDD option to the ccd item of the **Start Readout** command for higher imaging throughput.
- Changed calling convention to "stdcall" for compatibility with Visual C++ and Visual Basic.
- Changed calling struct member alignment to the default 8 bytes for compatibility.
- Made a single function callable from Windows 9X and Windows NT/2000/XP, using SBIGUDrvCommand() to replace ParDrvCommand() and ParDeviceCommandNH().
- Calls to **Open Device** and **Close Device** are now required.
- Added support for the USB and Ethernet based cameras and accessories.
- Removed the port information from the CC\_ESTABLISH\_LINK command and added it to the CC\_OPEN\_DEVICE command for compatibility with multiple hardware interfaces.
- Added the CC\_GET and CC\_SET\_DRIVER\_HANDLE commands for supporting multiple cameras per application.
- Added documentation for the **Get US Timer**, **Set IRQ**, **Get IRQ** commands.
- Added the **Get Line**, **Get Link Status**, **Start Readout** and **Get Error String** commands.
- Changed the minimum ST-7/8/9/10/1K exposure from 0.11 seconds to 0.12 seconds.

## Changes incorporated in Version 4.1

- Added readoutMode 9 with 9x9 binning to the ST-7/8/etc for faster focus mode throughput.

## Changes incorporated in Version 4.2

- Added support for the KAI2000 (ST-2K) and the TC-237 Tracking CCD
- Added requests 4 and 5 to the **Get CCD Info** command for use with the new frame transfer CCDs (TC237 Guider, KAI2000 Imager) when used with Start Exposure/End Exposure commands are bypassed as is typically done with AO exposures.
- Added the SBIGDriverChecker.exe Utility for updating the drivers to the current version.

## Changes incorporated in Version 4.21

- Added the **Set Driver Control** and **Get Driver Control** commands.
- Fixed a bug in the documentation where the ipAddress was ushort instead of ulong in the **Open Device** command.

## Changes incorporated in Version 4.22

- Added the **USB AD Control** command.
- Added the DCP\_CALL\_JOURNAL\_ENABLE and DCP\_IVTOH\_RATIO control parameters to the **Set Driver Control** command.
- Added the \Tools folder to the SBIG Universal Driver distribution with the EthSim.exe, SBIGUDRVJournalRx.exe, GetPortD.exe and SetClock.exe utility programs included.
- Released version 1.2 of the SBIGDriverChecker program that fixes a bug where some drivers could be installed in the wrong directory on some systems. Also removes those improperly installed drivers.
- Fix a bug in the documentation where the SBIGUnivDrvCommand() function was misnamed.
- Fix a bug in the documentation for the Activate Relay command regarding the bits in the command status.

## Changes incorporated in Version 4.23

- Fix many grammatical and typos with the much appreciated help of Andrew Mattingly.

## Changes incorporated in Version 4.24

- Documented the DCP\_IVTOH\_RATIO and added the DCP\_USB\_FIFO\_SIZE device control parameters which are Read or Set with the **Get/Set Driver Control** function.
- Extended the Exposure Range of the ST-2K with its Electronic Shutter from 0.12 seconds down to 0.01 seconds.

## Changes incorporated in Version 4.27

- added the **Query USB** command and the DEV\_USB1, DEV\_USB2, DEV\_USB3 and DEV\_USB4 enums to the Open Device command.

## Changes incorporated in Version 4.28

- Added the **Get Pentium Cycle Count** and **RW USB I2C** commands.

## Changes incorporated in Version 4.29/4.30

- Added the **CFW** command.
- Added support for the ST-L Large Format Camera. This mainly involved added support for the Remote Guiding Head CCD (CCD enum value 2)
- Added the **BitIO** command.

Changes incorporated in Version 4.30, 3<sup>rd</sup> Edition (Version 4.29 Build 16 of DLL)

- Made the CFW-L in the ST-L emulate a CFW-8 in the Pulse Out command.
- Added further support for the ST-L and CFW-L.
- Added support for the ST-402 camera and CFW-402 color filter wheel.
- Added the BITO\_FPGA\_WE sub-command to the Bit IO command.
- Added the DCP\_USB\_PIXEL\_DL\_ENABLE driver control parameter.
- Made one set of Driver Control Parameters for the driver vs. one per Handle.
- Added the CFWSEL\_AUTO select code for the **CFW** command.

Changes incorporated in Version 4.30, 4<sup>th</sup> Edition (Version 4.30 of DLL)

- Added the DCP\_HIGH\_THROUGHPUT driver control parameter.
- Allow 6 positions with the CFW-8 in the **CFW** command to support third party Color Filter Wheels that mimic the CFW-8 but have 6 filters.

## Changes incorporated in Version 4.35

- Documented the Update Clock function as inoperable.
- Add a new Section 2 describing the particulars of the Windows, Macintosh OSX and Linux versions of the library.
- Added the USB Loader Request to the Get Driver Info command for returning information about the USB Loader driver.
- Corrected the documentation regarding the ST-L's external shutter in the Remote Tracking Head via the **Miscellaneous Control** command. While the external shutter mimics the internal shutter for the Start Exposure command it does not with the Miscellaneous Control command. There you use separate control codes to affect the external shutter.

## Changes incorporated in Version 4.37

- Added support for millisecond exposures with Interline CCD cameras like the ST-2K, STL-11K and STL-4020. See the **Start Exposure** and **Get CCD Info** commands.

- For Single Shot Color cameras like the ST-2001XCM the name item in the **Get CCD Info** command contains the word “Color”. You can use this to detect whether to treat the data as Color or Monochrome for example.
- Added the DCP\_VDD\_OPTIMIZED parameter to the **Get/Set Device Control** commands.
- Fixed a bug where the DLL would not load under Windows 95.
- Added support for Trigger In capability with USB cameras. See the **Start Exposure** command.
- Added support for the very old CFW-6A to the **CFW** command.
- Added section 5.2 regarding Color CCDs.

#### Changes incorporated in Version 4.40

- Fixed a bug where CB\_CCD\_TYPE\_FRAME\_TRANSFER wasn't getting set for the STL-11K and STL-4020.
- Auto-select the USB driver (Standard or Alternate) based upon establishing a link.
- Added support for the KAI-2020 CCD version of the ST-2K
- Fixed a bug that would cause bright lines every 256 line in TDI mode.
- Added the DCP\_AUTO\_AD\_GAIN device control parameter.
- Added support for the STL-1301 model of the Large Format camera.
- Refer to the ST-402 by it's proper name instead of the development name (ST-F).

#### Changes incorporated in Version 4.42

- Added documentation for the ST-402, CFW-402 and CFW-10 products.

#### Changes incorporated in Version 4.42, 2<sup>nd</sup> Edition

- Refer to the CFW-402 instead of the development name CFW-F.

#### Changes incorporated in Version 4.42, 3<sup>rd</sup> Edition

- Expanded the documentation for the **CFW** Command including support for the CFWC\_OPEN\_DEVICE and CFWC\_CLOSE\_DEVICE sub-commands required for RS-232 control of the CFW-10.

#### Changes incorporated in Version 4.43

- Brought Linux and Macintosh drivers up to parity with the Windows version and released all as Version 4.43.
- No changes to API or documentation. Bug fixes in drivers only.

#### Changes incorporated in Version 4.44 Build 2

- Added the number of filter positions to the **cfwResult2** entry of the cfwResults struct for the CFWG\_FIRMWARE\_VERSION query of the CFWC\_GET\_INFO command.
- On cameras with a CFW-10 attached the Cfw Command will auto detect the CFW-10.
- Added support for the AOL, the AO for the STL to the **AO Tip Tilt** command.

#### Changes incorporated in Version 4.47

- Added support for the AOL (Large Format AO). The AOL is supported on the USB based ST-7/8/9/10/2K and STL series of cameras. The driver detects the type of AO attached to the camera in the Establish Link command and handles it accordingly. From the programming side the AOL acts just like an AO-7, but requires you to use the new AO Center command to recenter the AO.

- Added the Center AO command to force the AOL to re-center using the mechanical home sensors.
- Added the ability to disable the Auto Bias feature through the DriverChecker utility.
- Added support for the TDI mode on the interline Cameras like the ST-2000 and STL-11000.
- Added support for the external Remote Guide Head on the USB based ST-7/8/9/10/2K cameras. This includes setting flags in the capabilitiesBits field of the CCDInfoResults4 struct for Request 5 when an external Tracking Head is supported in hardware.

#### Changes incorporated in Version 4.52

- Added support for the KAI-11002 based STL-11K cameras.
- Added the ability to disable the vertical flush for KAI based cameras which gets used with ms level exposures to fix a smearing issue.
- Cleaned up and fixed a CFW-10 bug.
- Ported Mac versions to support the Intel based Macs.
- Fixed a bug in Serial based CFW-10 code.
- Fixed an intermittent Share error.
- Added ability to enable journaling through the Registry on Windows.
- Fixed a bug with the AO-L that would cause errors on large moves.
- Added support for the CFW-9.

#### Changes incorporated in Version 4.54

- Added support for the Biorad TDI mode of the ST-402.
- Added the **Motor Focus** command to support the Motor Focus accessories (section 5.7).
- Updated the help sections on supporting the new I<sup>2</sup>C based CFW and AO accessories.

#### Changes incorporated in Version 4.55

- Added preliminary support for the upcoming ST-H camera models.
- Added support for the low-cost AO (AO-8)

#### Changes incorporated in Version 4.57

- Added further support for the low-cost AO-8 (section 5.8).
- Added support for the ST-4000XCM, a small format 4-megapixel single shot color camera (section 5.9).

#### Changes incorporated in Version 4.60 Build 9

- Added the **Start Exposure 2** command (section 3.2.1) for use with the STX.
- Added the Advanced Results to the **Query Temperature Status** command (section 3.3.3).
- Added the **Query Ethernet** command for discovering Ethernet based SBIG Cameras (section 3.5.17).

#### Changes incorporated in Version 4.65 Build 3

- Added support for 64 bit versions of Windows to the DLL and Driver Checker.
- Added ability for the CFW-1602 to go to a special Dark Position with a Goto Position 4. After the move you must do a CFW Init before going to any other positions.
- Added several STX Production Code changes.
- Added support for the **SBIG BF2** format files for uploading firmware to the STX.

- Added **fanSpeed** and **trackingCCDSetpoint** to the **QueryTemperatureStatusResults2** struct for the **QueryTemperatureStatus** command. Also added the **TEMP\_STATUS\_ADVANCED2** **QUERY\_TEMP\_STATUS\_REQUEST** enum to retrieve this extra info.
- Added the **SetTemperatureRegulation2** command to allow setting the CCD Setpoint in degrees Celsius. See section 3.3.2.
- Added the **ReadOffset2** command. See section 3.5.8.
- Added several ST-8300 Production Changes.

#### Changes incorporated in Version 4.66

- Brought Mac OSX and PC Windows versions into parity.
- Added support for the STX-6303 model camera.
- Added developmental support for the Differential Guider accessory.
- Added support for Overscan pixels under Windows through a Registry switch.
- Additional production STX support.
- Added a bit to the Capabilities Bits in the Get CCD Info command that reports whether the camera has an internal memory based frame buffer like that found in the STX series of camera.

#### Changes incorporated in Version 4.67 Build 2

- Added the **ABORT\_DONT\_END** bit to the **ccd** field of the **End Exposure** command to allow completely terminating an STX image in progress.
- For an STX the **fanEnable** field can be used to manually control the fan speed of an STX as described in the **Miscellaneous Control** command notes.

#### Changes incorporated in Version 4.68 Build 4

- Added support for asynchronous (non-hanging) calls to the low level USB driver for the 32/64 bit Drivers which as accessible though the **Drv Options** button of the 64 Bit Driver Checker.
- Fix a bug under Windows 64 that could cause MaximDL to hang when establishing a link to a camera for the first time.
- Allowed a wider Overscan region on the ST-8300 cameras.
- Fixed a bug that would cause the FW\*-8300 CFW to only allow access to the first 5 filters instead of all eight.

#### Changes incorporated in Version 4.69 Build 23

- Added support for the ST-i model cameras.
- Deprecated the **Start Exposure** command. New and revised code should use the **Start Exposure 2** command instead for all cameras.

#### Changes incorporated in Version 4.70

- Bug fixes and internal testing software added.
- Simultaneous release under Windows and OS X.

#### Changes incorporated in Version 4.73B3

- Added preliminary support for the ST-8300B model camera.
- Improved performance with the ST-I model camera.
- Simultaneous release under Windows and OS X.

## Changes incorporated in Version 4.73B5

- Added Bit B6 to the capabilitiesBit field of the GetCCDInfo command which when set indicates the camera requires the StartExposure2 command instead of the older StartExposure command.

## Changes incorporated in Version 4.73B7

- Renamed the ST-8300B (the fast ST-8300) to the STF-8300.

## Changes incorporated in Version 4.74

- Added the DCP\_COLUMN\_REPAIR\_ENABLE, DCP\_WARM\_PIXEL\_REPAIR\_ENABLE, DCP\_WARM\_PIXEL\_REPAIR\_COUNT and DCP\_WARM\_PIXEL\_REPAIR\_COUNT Driver Control Parameters.
- Renamed a few cameras to be consistent with the web page.
- Fixed a bug with the ST-9 where the top few rows of long images could be black.
- Force the use of Asynchronous IO with the STF-8300 camera.

## Changes incorporated in Version 4.74B4

- Added the Column EEPROM command.
- Added the Fast Readout mode of the STF-8300 to the Start Exposure 2 command..

## Changes incorporated in Version 4.74B7

- Added Auto Filter support for the STX and STF-8300 cameras.
- February 2, 2012 - updated firmware for the STF-8300 and STX CCD cameras.

## Changes incorporated in Version 4.75B2

- Merged PC and OSX code into common code base.
- Added support for the FW-5 and the FW-7 (filter wheels) for the STX cameras.
- Added support for the FW-8 (filter wheel) for the STT cameras.

## Changes incorporated in Version 4.75B3

- Revision of the CFW documentation.

## Changes incorporated in Version 4.75B4

- Added support for CCD\_INFO\_REQUEST3 with enum value = 6, ie. GetCCDInfoResults6.

## Changes incorporated in Version 4.75B5

- Document improved.

## Changes incorporated in Version 4.75B6

- Added ST-I binning modes explanation to ReadoutLine function.

## Changes incorporated in Version 4.75B7

- Title page changed.

## Changes incorporated in Version 4.75B8

- Added Customer Options function.

## Changes incorporated in Version 4.75B9

- Added Debug Log function.

## Changes incorporated in Version 4.76B10

- Added support for up to eight USB and Ethernet based cameras.